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TMR (red) Tunel Cell Apoptosis Detection Kit

Product Description

Breakage of chromosomal DNA in apoptosis is a gradual phased process.

Chromosomal DNA is first degraded into large fragments of 50-300 kb by

endogenous nucleases, and then about 30% of chromosomal DNA is randomly

cleaved between nucleosome units under the action of Ca²⁺ and Mg²⁺-dependent

endonucleases to form 180-200 bp nucleosome DNA multimers. Therefore, in the late apoptosis stage, DNA is degraded into 180-200 bp fragments, exposing a large number of 3'-OH terminations on the broken genomic DNA. Terminal deoxynucleotidyl transferase (TdT) is a template-independent DNA polymerase that catalyzes deoxynucleotide binding to the 3'-OH end of a broken DNA molecule. Therefore, the TUNEL (TdT mediated dUTP Nick End Labeling) apoptosis detection kit can be used to detect the breakage of nuclear DNA in tissue cells during the advanced apoptosis process. The principle is that the 3'-OH end exposed at the time of genomic DNA breakage is incorporated with tetramethyl-rhodamine-dUTP (TMR-5-dUTP, TMR-5-dUTP) under the action of the TdT enzyme, which can be detected by fluorescence microscopy or flow cytometry (TMR excitation 551 nm, emission 575 nm). This kit has a wide range of applications, and is suitable for apoptosis detection of paraffin tissue sections, frozen tissue sections, cell crawling sections, cell smears, etc.

Storage and transportation

This kit is stored at -20°C, and TMR-5-dUTP Labeling Mix needs to be stored at -20°C in a dark place for 12 months.

Preparation Before the Experiment

1. PBS (Phosphate-Buffered Saline)
2. Fixative: 4% paraformaldehyde dissolved in PBS, pH 7.4
3. Permeabilization solution: 0.1%-0.5% Triton X-100

4. If nuclear staining is needed, prepare DAPI (2 $\mu\text{g}/\text{mL}$) yourself
5. If performing a positive control experiment, prepare DNase I yourself
6. If using a flow cytometer, prepare PI staining solution and RNase A (DNase free) yourself
7. Wear a lab coat and disposable gloves during operation.

Procedure:

1. Sample preparation

A. Paraffin-embedded tissue sections

1. Soak the paraffin tissue sections in environmentally friendly dewaxing transparent solution for 5-10 minutes at room temperature and repeat 3 times. Then soak in anhydrous ethanol for 5 minutes and repeat twice; Finally, soak in gradient ethanol (85% and 75%) and double distilled water once for 5 minutes each time.

2. Gently wash the sections with PBS and remove excess liquid around the sample. Use a grouping pen to draw a small circle with a distance of 2-3 mm from the tissue along the outer contour of the tissue to facilitate downstream permeability processing and balance marking operation. During the experiment, do not let the sample dry, and keep the sample moist in the wet box.

3. Preparation of Proteinase K working solution: According to the volume ratio of 1:9, PBS was used as a dilution to dilute Proteinase K (200 $\mu\text{g}/\text{mL}$) stock solution to make the final concentration 20 $\mu\text{g}/\text{mL}$.

4. Add 100 μ L of the above Proteinase K working solution to each sample, completely cover the tissue, and incubate at 37°C for 20 min;

(Note: Proteinase K treatment mainly helps the staining reagent to be permeable in the next steps of tissues and cells, and its incubation time will affect the subsequent labeling efficiency if it is too long or too short, so the incubation time of Proteinase K can be optimized for better results)

5. Wash the samples with PBS solution infiltrate 3 times for 5 min each time (Proteinase K needs to be washed clean, otherwise it will interfere with subsequent labeling reactions). The treated samples are placed in a wet box to keep the samples moist;

6. (Optional step) Remove the excess liquid from the sample, add an appropriate amount of membrane rupture droplets to the tissue, fully infiltrate the tissue, and treat at room temperature for 20 minutes. After the film breaking treatment was completed, the samples were washed with PBS solution 3 times for 5 minutes each time. The treated sample is placed in a wet box to keep the sample moist.

B. Organizing Frozen Sections

1. Immerse the slides in 4% paraformaldehyde solution (dissolved in PBS) for fixation and incubate at room temperature for 10-15 minutes;

2. Remove the slides from the fixative and air dry them naturally in a fume hood;

3. Rinse the slides in pure water or PBS to remove any residual fixative on the slides;

4. Use a tissue-marking pen to draw a small circle around the tissue outline, leaving

a 2-3 mm gap from the tissue, to facilitate downstream permeability treatment and balanced labeling operations; during the experiment, do not let the samples dry, and keep treated samples in a humid box to maintain moisture;

5. Prepare Proteinase K working solution: dilute the Proteinase K stock solution (200 $\mu\text{g}/\text{mL}$) with PBS at a ratio of 1:9 to reach a final concentration of 20 $\mu\text{g}/\text{mL}$;

6. Add 100 μL of the Proteinase K working solution to each sample, ensuring complete coverage, and incubate at room temperature for 10 minutes;

(Note: Proteinase K treatment mainly helps improve the permeability of tissues and cells for subsequent staining reagents. Incubation time that is too long or too short may affect the efficiency of downstream labeling. If optimal results are not achieved, the Proteinase K incubation time may need to be optimized.)

7. Rinse the samples 2-3 times with PBS to remove excess liquid (Proteinase K must be thoroughly washed off, otherwise it will interfere with subsequent labeling reactions). Keep the treated samples in a humid box to maintain moisture;

8. (Optional) Add an appropriate amount of permeabilization solution to the tissue, fully infiltrate the tissue, and treat at room temperature for 20 minutes. After permeabilization, rinse the samples with PBS to remove excess liquid, and keep the treated samples in a humid box to maintain moisture.

C. Cell Spreading on Coverslips

1. Culture adherent cells in confocal dishes. After apoptosis induction treatment, gently rinse the cells twice with PBS;

2. Add an appropriate amount of 4% paraformaldehyde solution (dissolved in PBS) to each small dish for fixation and incubate at room temperature for 20 minutes;
3. Remove the fixative and wash three times with PBS, 5 minutes each time;
4. Immerse each sample in permeabilization solution and incubate at room temperature for 5 minutes for permeabilization (Note: It is recommended to digest with 2-20 $\mu\text{g/mL}$ Proteinase K working solution at 37°C for about 10 minutes, adjusting as needed based on cell condition. If cells are prone to detachment, using permeabilization solution is recommended);
5. Immerse the samples in an open beaker containing PBS solution and rinse 2-3 times;
6. Gently remove excess liquid and carefully blot the surroundings of the samples on the slides with filter paper. Place the treated samples in a humidified box to keep them moist.

D. Cell Smear

1. Resuspend the cells in PBS at a concentration of approximately 2×10^7 cells/mL. Pipette 50-100 μL of the cell suspension onto an anti-fall glass slide, and gently spread the cell suspension using a clean coverslip.
2. Immerse the slide in a staining jar containing 4% freshly prepared paraformaldehyde in PBS to fix the cells, and leave it at 4°C for 25 minutes.
3. Immerse the slide in PBS and let it wash at room temperature for 5 minutes, then repeat once.

4. Gently remove the excess liquid and carefully absorb any remaining liquid around the sample on the slide with filter paper. Draw a small circle around the cell perimeter with a marker to facilitate subsequent permeability treatment and labeling procedures. During the experiment, do not let the sample dry.

5. Immerse each sample in permeabilization solution and incubate at room temperature for 5 minutes for permeabilization (Note: It is recommended to digest with Proteinase K working solution at 2-20 $\mu\text{g}/\text{mL}$, treating at 37°C for about 10 minutes, adjusting according to the cell condition. If cells tend to detach, it is recommended to use permeabilization solution instead).

6. Immerse and wash the samples 2-3 times in an open beaker containing PBS solution.

7. Gently remove the excess liquid and carefully absorb any remaining liquid around the sample on the slide with filter paper. Place the treated samples in a humidified chamber to keep them moist.

2. DNase I Treatment Positive Control Experiment (Optional Step)

After sample permeabilization, treat the samples with DNase I to prepare a positive control.

1. Add 100 μL of $1\times$ DNase I Buffer (prepared by mixing 10 μL of $10\times$ DNase I Buffer with 90 μL of deionized water) to the permeabilized samples and incubate at room temperature for 5 minutes;

2. Gently remove the excess liquid, add 100 μL of the working solution containing

DNase I (20 U/mL), and incubate at room temperature for 10 minutes;

3. Gently remove the excess liquid and thoroughly wash the slides 3–4 times in a staining jar containing PBS.

3. Labeling and Detection

1. Equilibration: Add 100 μ L of Equilibration Buffer to each sample to completely cover the area to be tested, and incubate at room temperature for 10 minutes.

2. Labeling solution preparation: Thaw TMR-5-dUTP Labeling Mix and Equilibration Buffer on ice, and mix them according to the ratio Recombinant TdT enzyme:TMR-5-dUTP Labeling Mix:Equilibration Buffer = 2 μ L:10 μ L:100 μ L (2:10:100) to make enough TdT incubation buffer for all experiments. The volume of reagents used in specific experiments can be appropriately adjusted proportionally based on the size of the slides.

3. Negative control system: Prepare a control TdT incubation buffer without Recombinant TdT enzyme, replacing it with ddH₂O.

4. Labeling: Remove as much of the equilibrated Equilibration Buffer as possible, then add 112 μ L of TdT incubation buffer to each tissue sample and incubate at 37°C for 1 hour. Ensure the slides do not dry out and keep them away from light.

5. Immediately wash the tissue samples with PBS three times, 5 minutes each wash.

6. Gently wipe off PBS around the samples with filter paper.

7. Nuclear staining: Place the samples in a staining jar, immerse the slides in a jar containing DAPI solution (freshly prepared and diluted with PBS) in the dark, and

let stand at room temperature for 8 minutes.

8. Mounting: After staining, wash the tissue samples with PBS three times, 5 minutes each time, gently remove excess liquid, and add anti-fade mounting medium for mounting.

9. Microscopy: Analyze the samples immediately under a fluorescence microscope, keeping the slides away from light. DAPI stains both apoptotic and non-apoptotic cells blue, while the red fluorescence from TMR-5-dUTP incorporation is only localized in the nuclei of apoptotic cells.

4. Detection of suspended cells by flow cytometry

1. Wash the cells to be tested twice with PBS, centrifuge at 4 °C (500 g) and resuspend in 500 μL PBS;

2. Fixation: Add 5 mL of 1% paraformaldehyde solution prepared with PBS to the sample, fix the cells, and place on ice for 20 minutes;

3. Cells were centrifuged at 300 g at 4 °C for 10 min, desupernatant and resuspended twice with 5 mL PBS, and finally resuspended with 500 μL PBS.

4. Permeabilization: Add 5 mL of pre-chilled 70% ethanol on ice to the sample, incubate at -20°C for 4 h, permeate the cells;

(Note: Cells can also be permeated by membrane breaking solution at room temperature for 5 minutes)

5. Cells were centrifuged with 300 g for 10 min and resuspended with 5 mL PBS, and re-centrifuged with 1 mL PBS.

6. Equilibration: Transfer approximately 2×10^6 cells to 1.5 mL microcentrifuge tubes, centrifuge 300 g for 10 min, remove the supernatant, and resuspend with 100 μ L Equilibration Buffer and incubate for 5 min at room temperature.
7. Labeling solution preparation: Thaw TMR-5-dUTP Labeling Mix and Equilibration Buffer on ice and follow Recombinant TdT enzyme: TMR-5-dUTP Labeling Mix: Equilibration Buffer=2 μ L: 10 μ L: 100 μ L (2:10:100) ratio mix enough TdT incubation buffer for all experiments and optional positive control reactions;
8. Labeling: Cells were centrifuged with 300 g for 10 min, desupernatant and resuspended in 112 μ L TdT incubation buffer and incubated at 37 $^{\circ}$ C for 1 h away from light. Gently resuspend cells with a micropipette every 15 min;
9. After the reaction is completed, add 1 mL of 20 mM EDTA to terminate the reaction and gently mix well with a micropipette.
10. Centrifuge 300 g for 10 min, remove the supernatant and resuspend the pellet in 1 mL of membrane breaking solution containing 5 mg/mL BSA, repeat the wash twice;
11. Nuclear staining: 300 g centrifuge for 10 min, desupernatant and resuspend the cell pellet in 0.5 mL DAPI solution containing 250 μ g DNase-free RNase A, incubate cells at room temperature in the dark for 30 min;
12. On-machine detection: Flow cytometry analyzes cells, DAPI can stain both apoptotic and non-apoptotic nuclei in blue, and only the apoptotic nuclei have red fluorescence localized by TMR-5-dUTP incorporation.