

## Anti-RTKN2 antibody

**Cat. No.** ml125757  
**Package** 25 µl/100 µl/200 µl  
**Storage** -20°C, pH7.4 PBS, 0.05% NaN<sub>3</sub>, 40% Glycerol

### Product overview

**Description** Anti-RTKN2 rabbit polyclonal antibody  
**Applications** ELISA, IHC  
**Immunogen** Fusion protein of human RTKN2  
**Reactivity** Human  
**Content** 1.2 mg/ml  
**Host species** Rabbit  
**Ig class** Immunogen-specific rabbit IgG  
**Purification** Antigen affinity purification

### Target information

**Symbol** RTKN2  
**Full name** rhotekin 2  
**Synonyms** PLEKHK1; bA531F24.1  
**Swissprot** Q8IZC4

### Target Background

May play an important role in lymphopoiesis.

订购热线: 4008-898-798


#### Applications

##### Immunohistochemistry

Predicted cell location: Cytoplasm

Positive control: Human liver cancer

Recommended dilution: 50-300



The image on the left is immunohistochemistry of paraffin-embedded Human liver cancer tissue using ml125757(RTKN2 Antibody) at dilution 1/70, on the right is treated with fusion protein. (Original magnification:  $\times 200$ )

##### ELISA

Recommended dilution: 5000-10000

联系电话: 4008-898-798, 021-61725725

联系QQ: 2881505695, 2881505696

邮箱: [mlbio\\_cn@yeah.net](mailto:mlbio_cn@yeah.net)

网址: [www.mlbio.cn](http://www.mlbio.cn)