

订购热线: 4008-898-798

## Anti-TEX11 antibody

**Cat. No.** ml225527

**Package** 25 μl/100 μl/200 μl

**Storage** -20°C, pH7.4 PBS, 0.05% NaN3, 40% Glycerol

**Product overview** 

**Description** Anti-TEX11 rabbit polyclonal antibody

**Applications** ELISA, IHC

**Immunogen** Fusion protein of human TEX11

ReactivityHumanContent0.54 mg/mlHost speciesRabbit

Ig classImmunogen-specific rabbit IgGPurificationAntigen affinity purification

**Target information** 

Symbol TEX11

Full nametestis expressed 11SynonymsTGC1; TSGA3; SPGFX2

Swissprot Q8IYF3

## **Target Background**

This gene is X-linked and is expressed in only male germ cells. Two alternatively spliced transcript variants encoding distinct isoforms have been found for this gene.



订购热线: 4008-898-798

## **Applications**

## **Immunohistochemistry**

Predicted cell location: Cytoplasm and Nucleus

Positive control: Human lung cancer Recommended dilution: 30-150





The image on the left is immunohistochemistry of paraffin-embedded Human lung cancer tissue using ml225527(TEX11 Antibody) at dilution 1/30, on the right is treated with fusion protein. (Original magnification: ×200)

**ELISA** 

Recommended dilution: 5000-10000

联系电话: 4008-898-798, 021-61725725


联系QQ: 2881505695, 2881505696

邮箱: mlbio\_cn@yeah.net 网址: www.mlbio.cn

Predicted cell location: Cytoplasm and Nucleus Positive control: Human esophagus cancer

Recommended dilution: 30-150





The image on the left is immunohistochemistry of paraffin-embedded Human esophagus cancer tissue using ml225527(TEX11 Antibody) at dilution 1/30, on the right is treated with fusion protein. (Original magnification: ×200)