

订购热线: 4008-898-798

Anti-SLC41A2 antibody

Cat. No. ml260743

Package 25 μl/100 μl/200 μl

Storage -20°C, pH7.4 PBS, 0.05% NaN3, 40% Glycerol

Product overview

Description Anti-SLC41A2 rabbit polyclonal antibody

Applications ELISA, IHC

Immunogen Synthetic peptide of human SLC41A2

ReactivityHumanContent0.5 mg/mlHost speciesRabbit

Ig classImmunogen-specific rabbit IgGPurificationAntigen affinity purification

Target information

Symbol SLC41A2

Full name Solute carrier family 41, member 2

Synonyms SLC41A1-L1
Swissprot Q96JW4

Target Background

This gene belongs to the SLC41A transporter family. It acts as a plasma-membrane magnesium transporter. SLC41A2 functions as a plasma-membrane Mg2+ transporter in vertebrate cells.





订购热线: 4008-898-798

Applications

Immunohistochemistry

Predicted cell location: Cytoplasm Positive control: Human gasrtic cancer Recommended dilution: 25-100





The image on the left is immunohistochemistry of paraffin-embedded Human gasrtic cancer tissue using ml260743(SLC41A2 Antibody) at dilution 1/30, on the right is treated with synthetic peptide. (Original magnification: ×200)

ELISA

Recommended dilution: 1000-5000

联系电话: 4008-898-798, 021-61725725

联系QQ: 2881505695, 2881505696

邮箱: mlbio_cn@yeah.net 网址: www.mlbio.cn