

订购热线: 4008-898-798

Anti-SLC34A2 antibody

Cat. No. ml260744

Package 25 μ l/100 μ l/200 μ l

Storage -20°C, pH7.4 PBS, 0.05% NaN3, 40% Glycerol

Product overview

Description Anti-SLC34A2 rabbit polyclonal antibody

Applications ELISA, IHC

Immunogen Synthetic peptide of human SLC34A2

ReactivityHumanContent0.6 mg/mlHost speciesRabbit

Ig classImmunogen-specific rabbit IgGPurificationAntigen affinity purification

Target information

Symbol SLC34A2

Full name Solute carrier family 34 (sodium phosphate), member 2

Synonyms NPTIIb, NAPI-3B, NAPI-IIb

Swissprot O95436

Target Background

The protein encoded by this gene is a pH-sensitive sodium-dependent phosphate transporter. Phosphate uptake is increased at lower pH. Defects in this gene are a cause of pulmonary alveolar microlithiasis. Three transcript variants encoding two different isoforms have been found for this gene.




订购热线: 4008-898-798

Applications

Immunohistochemistry


Predicted cell location: Cytoplasm Positive control: Human gasrtic cancer Recommended dilution: 50-200






The image on the left is immunohistochemistry of paraffin-embedded Human gasrtic cancer tissue using ml260744(SLC34A2 Antibody) at dilution 1/50, on the right is treated with synthetic peptide. (Original magnification: ×200)

Predicted cell location: Cytoplasm Positive control: Human ovarian cancer Recommended dilution: 50-200





The image on the left is immunohistochemistry of paraffin-embedded Human ovarian cancer tissue using ml260744(SLC34A2 Antibody) at dilution 1/50, on the right is treated with synthetic peptide. (Original magnification: ×200)

ELISA

Recommended dilution: 2000-10000

联系电话: 4008-898-798, 021-61725725

联系QQ: 2881505695, 2881505696

邮箱: mlbio_cn@yeah.net 网址: www.mlbio.cn