

Anti-ABCB5 antibody

Cat. No.	ml120189
Package	25 µl/100 µl/200 µl
Storage	-20°C, pH7.4 PBS, 0.05% NaN ₃ , 40% Glycerol

Product overview

Description	Anti-ABCB5 rabbit polyclonal antibody
Applications	ELISA, IHC
Immunogen	Fusion protein of human ABCB5
Reactivity	Human, Mouse
Content	0.4 mg/ml
Host species	Rabbit
Ig class	Immunogen-specific rabbit IgG
Purification	Antigen affinity purification

Target information

Symbol	ABCB5
Full name	ATP-binding cassette sub-family B member 5
Synonyms	ABCB5beta, EST422562, ABCB5alpha
Swissprot	Q2M3G0

Target Background

ABCB5 belongs to the ATP-binding cassette (ABC) transporter superfamily of integral membrane proteins. These proteins participate in ATP-dependent transmembrane transport of structurally diverse molecules ranging from small ions, sugars, and peptides to more complex organic molecules.

订购热线: 4008-898-798


Applications

Immunohistochemistry

Predicted cell location: Cell membrane, Cytoplasm

Positive control: Human cervical cancer

Recommended dilution: 10-50




The image is immunohistochemistry of paraffin-embedded Human cervical cancer tissue using ml120189(ABCB5 Antibody) at dilution 1/10. (Original magnification: ×200)

Predicted cell location: Cell membrane, Cytoplasm

Positive control: Human renal cancer

Recommended dilution: 10-50



The image is immunohistochemistry of paraffin-embedded Human renal cancer tissue using ml120189(ABCB5 Antibody) at dilution 1/10. (Original magnification: ×200)

ELISA

Recommended dilution: 500-2000

联系电话: 4008-898-798, 021-61725725

联系QQ: 2881505695, 2881505696

邮箱: mlbio_cn@yeah.net

网址: www.mlbio.cn