

订购热线: 4008-898-798

Anti-MDK antibody

Cat. No. ml120984

Package 25 μl/100 μl/200 μl

Storage -20°C, pH7.4 PBS, 0.05% NaN3, 40% Glycerol

Product overview

Description Anti-MDK rabbit polyclonal antibody

Applications ELISA, IHC

Immunogen Fusion protein of human MDK

ReactivityHumanContent0.2 mg/mlHost speciesRabbit

Ig classImmunogen-specific rabbit IgGPurificationAntigen affinity purification

Target information

Symbol MDK

Full name midkine (neurite growth-promoting factor 2)

Synonyms MK, ARAP, NEGF2

Swissprot P21741

Target Background

This gene encodes a member of a small family of secreted growth factors that binds heparin and responds to retinoic acid. The encoded protein promotes cell growth, migration, and angiogenesis, in particular during tumorigenesis. This gene has been targeted as a therapeutic for a variety of different disorders. Alternatively spliced transcript variants encoding multiple isoforms have been observed.




订购热线: 4008-898-798

Applications

Immunohistochemistry


Predicted cell location: Cytoplasm Positive control: Human colon cancer Recommended dilution: 50-200






The image on the left is immunohistochemistry of paraffin-embedded Human colon cancer tissue using ml120984(MDK Antibody) at dilution 1/60, on the right is treated with fusion protein. (Original magnification: ×200)

Predicted cell location: Cytoplasm Positive control: Human lung cancer Recommended dilution: 50-200





The image on the left is immunohistochemistry of paraffin-embedded Human lung cancer tissue using ml120984(MDK Antibody) at dilution 1/60, on the right is treated with fusion protein. (Original magnification: ×200)

ELISA

Recommended dilution: 2000-5000

联系电话: 4008-898-798, 021-61725725

联系QQ: 2881505695, 2881505696

邮箱: mlbio_cn@yeah.net 网址: www.mlbio.cn