

Anti-ARNTL antibody

Cat. No.	ml121991
Package	25 µl/100 µl/200 µl
Storage	-20°C, pH7.4 PBS, 0.05% NaN ₃ , 40% Glycerol

Product overview

Description	Anti-ARNTL rabbit polyclonal antibody
Applications	ELISA, WB
Immunogen	Fusion protein of human ARNTL
Reactivity	Human, Mouse, Rat
Content	0.9 mg/ml
Host species	Rabbit
Ig class	Immunogen-specific rabbit IgG
Purification	Antigen affinity purification

Target information

Symbol	ARNTL
Full name	aryl hydrocarbon receptor nuclear translocator-like
Synonyms	TIC; JAP3; MOP3; BMAL1; PASD3; BMAL1c; bHLHe5
Swissprot	O00327

Target Background

The protein encoded by this gene is a basic helix-loop-helix protein that forms a heterodimer with CLOCK. This complex binds an E-box upstream of the PER1 gene, activating this gene and possibly other circadian rhythm-associated genes. Three transcript variants encoding two different isoforms have been found for this gene.

订购热线: 4008-898-798

Applications

Western blotting

Predicted band size:69 kDa

Positive control:Human fetal liver and mouse stomach tissue, hela and hepG2 cells

Recommended dilution: 500-2000

Gel: 6%SDS-PAGE


Lysate: 40 μ g

Lane 1-4: Human fetal liver tissue, mouse stomach tissue,
hela cells, hepG2 cells

Primary antibody: ml121991(ARNTL Antibody) at dilution
1/450

Secondary antibody: Goat anti rabbit IgG at 1/8000 dilution

Exposure time: 20 seconds



ELISA

Recommended dilution:2000-5000

联系电话: 4008-898-798, 021-61725725

联系QQ: 2881505695, 2881505696

邮箱: mlbio_cn@yeah.net

网址: www.mlbio.cn