

订购热线: 4008-898-798

Anti-CTNNBIP1 antibody

Cat. No. ml122164

Package 25 μl/100 μl/200 μl

Storage -20°C, pH7.4 PBS, 0.05% NaN3, 40% Glycerol

Product overview

Description Anti-CTNNBIP1 rabbit polyclonal antibody

Applications ELISA, IHC

Immunogen Fusion protein of human CTNNBIP1

ReactivityHuman, MouseContent0.4 mg/mlHost speciesRabbit

Ig classImmunogen-specific rabbit IgGPurificationAntigen affinity purification

Target information

Symbol CTNNBIP1

Full name catenin, beta interacting protein 1

Synonyms ICAT
Swissprot Q9NSA3

Target Background

The protein encoded by this gene binds CTNNB1 and prevents interaction between CTNNB1 and TCF family members. The encoded protein is a negative regulator of the Wnt signaling pathway. Two transcript variants encoding the same protein have been found for this gene.





订购热线: 4008-898-798

Applications

Immunohistochemistry


Predicted cell location: Cytoplasm Positive control: Human breast cancer Recommended dilution: 50-200






The image on the left is immunohistochemistry of paraffin-embedded Human breast cancer tissue using ml122164(CTNNBIP1 Antibody) at dilution 1/30, on the right is treated with fusion protein. (Original magnification: ×200)

Predicted cell location: Cytoplasm Positive control: Human thyroid cancer Recommended dilution: 50-200





The image on the left is immunohistochemistry of paraffin-embedded Human thyroid cancer tissue using ml122164(CTNNBIP1 Antibody) at dilution 1/30, on the right is treated with fusion protein. (Original magnification: ×200)

ELISA

Recommended dilution: 2000-5000

联系电话: 4008-898-798, 021-61725725

联系QQ: 2881505695, 2881505696

邮箱: mlbio_cn@yeah.net 网址: www.mlbio.cn