

Anti-UTS2R antibody

Cat. No.	ml160449
Package	25 µl/100 µl/200 µl
Storage	-20°C, pH7.4 PBS, 0.05% NaN ₃ , 40% Glycerol

Product overview

Description	Anti-UTS2R rabbit polyclonal antibody
Applications	ELISA, WB
Immunogen	Synthetic peptide of human UTS2R
Reactivity	Human
Content	0.2 mg/ml
Host species	Rabbit
Ig class	Immunogen-specific rabbit IgG
Purification	Antigen affinity purification

Target information

Symbol	UTS2R
Full name	urotensin 2 receptor
Synonyms	UTR, UTR2, GPR14, UR-2-R
Swissprot	Q9UKP6

Target Background

This gene codes high affinity receptor for urotensin-2 and urotensin-2B. The activity of this receptor is mediated by a G-protein that activate a phosphatidylinositol-calcium second messenger system

订购热线: 4008-898-798

Applications

Western blotting

Predicted band size:42 kDa

Positive control:Human fetal muscle tissue

Recommended dilution: 500-2000

Gel: 10%SDS-PAGE


Lysate: 60 µg

Lane: Human fetal muscle tissue

Primary antibody: ml160449(UTS2R Antibody) at dilution 1/400

Secondary antibody: Goat anti rabbit IgG at 1/8000 dilution

Exposure time: 30 seconds



ELISA

Recommended dilution:1000-5000

联系电话: 4008-898-798, 021-61725725

联系QQ: 2881505695, 2881505696

邮箱: mlbio_cn@yeah.net

网址: www.mlbio.cn