

Anti-SOD2 antibody

Cat. No.	ml221247
Package	25 µl/100 µl/200 µl
Storage	-20°C, pH7.4 PBS, 0.05% NaN ₃ , 40% Glycerol

Product overview

Description	Anti-SOD2 rabbit polyclonal antibody
Applications	ELISA, WB
Immunogen	Fusion protein of human SOD2
Reactivity	Human, Mouse, Rat
Content	0.2 mg/ml
Host species	Rabbit
Ig class	Immunogen-specific rabbit IgG
Purification	Antigen affinity purification

Target information

Symbol	SOD2
Full name	Superoxide dismutase 2, mitochondrial
Synonyms	IPOB; MNSOD; MVCD6
Swissprot	P04179

Target Background

This gene is a member of the iron/manganese superoxide dismutase family. It encodes a mitochondrial protein that forms a homotetramer and binds one manganese ion per subunit. This protein binds to the superoxide byproducts of oxidative phosphorylation and converts them to hydrogen peroxide and diatomic oxygen. Mutations in this gene have been associated with idiopathic cardiomyopathy (IDC), premature aging, sporadic motor neuron disease, and cancer. Alternate transcriptional splice variants, encoding different isoforms, have been characterized.

订购热线: 4008-898-798

Applications

Western blotting

Predicted band size:25 kDa

Positive control:Mouse liver tissue

Recommended dilution: 500-2000

Gel: 12%SDS-PAGE


Lysate: 40 μ g

Lane: Mouse liver tissue

Primary antibody: ml221247(SOD2 Antibody) at dilution 1/650

Secondary antibody: Goat anti rabbit IgG at 1/8000 dilution

Exposure time: 30 seconds



ELISA

Recommended dilution:2000-5000

联系电话: 4008-898-798, 021-61725725

联系QQ: 2881505695, 2881505696

邮箱: mlbio_cn@yeah.net

网址: www.mlbio.cn