

订购热线: 4008-898-798

Anti-TYR antibody

Cat. No.	ml221392
Package	25 μl/100 μl/200 μl
Storage	-20°C, pH7.4 PBS, 0.05% NaN3, 40% Glycerol
Product overview	
Description	Anti-TYR rabbit polyclonal antibody
Applications	ELISA, WB, IHC
Immunogen	Fusion protein of human TYR
Reactivity	Human, Mouse
Content	0.2 mg/ml
Host species	Rabbit
lg class	Immunogen-specific rabbit IgG
Purification	Antigen affinity purification
	Antigen affinity purification
Target information	oroc
Symbol	TYR
Full name	Tyrosinase
(,00



Synonyms CMM8; OCA1A; OCAIA; SHEP3

Swissprot P14679

Target Background

The enzyme encoded by this gene catalyzes the first 2 steps, and at least 1 subsequent step, in the conversion of tyrosine to melanin. The enzyme has both tyrosine hydroxylase and dopa oxidase catalytic activities, and requires copper for function. Mutations in this gene result in oculocutaneous albinism, and nonpathologic polymorphisms result in skin pigmentation variation. The human genome contains a pseudogene similar to the 3' half of this gene.




Applications

Immunohistochemistry

Predicted cell location: Cytoplasm

Positive control: Human thyroid cancer

Recommended dilution: 25-100



The image on the left is immunohistochemistry of paraffin-embedded Human thyroid cancer tissue using ml221392(TYR Antibody) at dilution 1/30, on the right is treated with fusion protein. (Original magnification: ×200)

Good elisa Nil producer

Western blotting

Predicted band size:60 kDa


Positive control:A375, Hela and Jurkat cells


Recommended dilution: 500-2000

Predicted cell location: Cytoplasm

Positive control: Human lung cancer

Recommended dilution: 25-100





The image on the left is immunohistochemistry of paraffin-embedded Human lung cancer tissue using ml221392(TYR Antibody) at dilution 1/30, on the right is treated with fusion protein. (Original magnification: ×200)



Gel: 10%SDS-PAGE	kDa 1 2 3 130 —
Lysate: 50 µg	130 — 100 —
Lane 1-3: A375 cells, Hela cells, Jurkat cells	70 — 55 —
Primary antibody: ml221392(TYR Antibody) at dilution 1/400	40 —
Secondary antibody: Goat anti rabbit IgG at 1/8000 dilution	35 —
Exposure time: 3 minutes	25 —

ELISA

Recommended dilution: 2000-5000

联系电话: 4008-898-798,021-61725725

联系QQ: 2881505695, 2881505696

邮箱: mlbio_cn@yeah.net

网址: www.mlbio.cn