

Anti-CERKL antibody

Cat. No.	ml261556
Package	25 µl/100 µl/200 µl
Storage	-20°C, pH7.4 PBS, 0.05% NaN ₃ , 40% Glycerol

Product overview

Description	Anti-CERKL rabbit polyclonal antibody
Applications	ELISA, WB
Immunogen	Synthetic peptide of human CERKL
Reactivity	Human
Content	0.4 mg/ml
Host species	Rabbit
Ig class	Immunogen-specific rabbit IgG
Purification	Antigen affinity purification

Target information

Symbol	CERKL
Full name	ceramide kinase-like
Synonyms	RP26
Swissprot	Q49MI3

Target Background

This gene encodes a protein with ceramide kinase-like domains, however, the protein does not phosphorylate ceramide and its target substrate is currently unknown. This protein may be a negative regulator of apoptosis in photoreceptor cells. Mutations in this gene cause a form of retinitis pigmentosa characterized by autosomal recessive cone and rod dystrophy (arCRD).

订购热线: 4008-898-798

Applications

Western blotting

Predicted band size:63 kDa

Positive control:Hela cells

Recommended dilution: 500-2000

Gel: 6%SDS-PAGE


Lysate: 40 µg

Lane: Hela cells

Primary antibody: ml261556(CERKL Antibody) at dilution 1/1100

Secondary antibody: Goat anti rabbit IgG at 1/8000 dilution

Exposure time: 1 minute



ELISA

Recommended dilution:2000-5000

联系电话: 4008-898-798, 021-61725725

联系QQ: 2881505695, 2881505696

邮箱: mlbio_cn@yeah.net

网址: www.mlbio.cn