

Anti-NANOG antibody

Cat. No.	ml262945
Package	25 µl/100 µl/200 µl
Storage	-20°C, pH7.4 PBS, 0.05% NaN ₃ , 40% Glycerol

Product overview

Description	Anti-NANOG rabbit polyclonal antibody
Applications	ELISA, WB
Immunogen	Synthetic peptide of human NANOG
Reactivity	Human, Mouse
Content	0.1 mg/ml
Host species	Rabbit
Ig class	Immunogen-specific rabbit IgG
Purification	Antigen affinity purification

Target information

Symbol	NANOG
Full name	Nanog homeobox
Synonyms	
Swissprot	Q9H9S0

Target Background

The protein encoded by this gene is a DNA binding homeobox transcription factor involved in embryonic stem (ES) cell proliferation, renewal, and pluripotency. The encoded protein can block ES cell differentiation and can also autorepress its own expression in differentiating cells. Two transcript variants encoding different isoforms have been found for this gene.

订购热线: 4008-898-798

Applications

Western blotting

Predicted band size:35 kDa

Positive control:Mouse brain tissue

Recommended dilution: 200-1000

Gel: 12%SDS-PAGE


Lysate: 40 µg

Lane: Mouse brain tissue

Primary antibody: ml262945(NANOG Antibody) at dilution 1/200

Secondary antibody: Goat anti rabbit IgG at 1/8000 dilution

Exposure time: 5 minutes



ELISA

Recommended dilution:5000-10000

联系电话: 4008-898-798, 021-61725725

联系QQ: 2881505695, 2881505696

邮箱: mlbio_cn@yeah.net

网址: www.mlbio.cn