

订购热线: 4008-898-798

Anti-STEAP2 antibody

Cat. No. ml262955

Package 25 μ l/100 μ l/200 μ l

Storage -20°C, pH7.4 PBS, 0.05% NaN3, 40% Glycerol

Product overview

Description Anti-STEAP2 rabbit polyclonal antibody

Applications ELISA, WB

Immunogen Synthetic peptide of human STEAP2

ReactivityHuman, MouseContent0.9 mg/mlHost speciesRabbit

Ig classImmunogen-specific rabbit IgGPurificationAntigen affinity purification

Target information

Symbol STEAP2

Full nameSTEAP family member 2, metalloreductaseSynonymsSTMP; IPCA1; PUMPCn; STAMP1; PCANAP1

Swissprot Q8NFT2

Target Background

This gene is a member of the STEAP family and encodes a multi-pass membrane protein that localizes to the Golgi complex, the plasma membrane, and the vesicular tubular structures in the cytosol. A highly similar protein in mouse has both ferrireductase and cupric reductase activity, and stimulates the cellular uptake of both iron and copper in vitro. Increased transcriptional expression of the human gene is associated with prostate cancer progression. Alternate transcriptional splice variants, encoding different isoforms, have been characterized.



订购热线: 4008-898-798

Applications

Western blotting

Predicted band size:56 kDa

Positive control:Human prostate tissue Recommended dilution: 500-2000

Gel: 8%SDS-PAGE

Lysate: 40 µg


Lane: Human prostate tissue

Primary antibody: ml262955(STEAP2 Antibody) at dilution

1/400

Secondary antibody: Goat anti rabbit IgG at 1/8000 dilution

Exposure time: 20 seconds



ELISA

Recommended dilution:5000-10000

联系电话: 4008-898-798, 021-61725725

联系QQ: 2881505695, 2881505696

邮箱: mlbio_cn@yeah.net

网址: www.mlbio.cn