

# Anti-SLC17A8 antibody

Cat. No.	ml263210
Package	25 μl/100 μl/200 μl
Storage	-20°C, pH7.4 PBS, 0.05% NaN3, 40% Glycerol

Product overview	
Description	Anti-SLC17A8 rabbit polyclonal antibody
Applications	ELISA, WB
Immunogen	Synthetic peptide of human SLC17A8
Reactivity	Human
Content	0.6 mg/ml
Host species	Rabbit
lg class	Immunogen-specific rabbit IgG
Purification	Antigen affinity purification
Target information	
Symbol	SLC17A8

30000

ranget internation	
Symbol	SLC17A8
Full name	solute carrier family 17 (vesicular glutamate transporter), member 8
Synonyms	DFNA25; VGLUT3
Swissprot	Q8NDX2
Target Background	odulo odulo

### **Target Background**


This gene encodes a vesicular glutamate transporter. The encoded protein transports the neurotransmitter glutamate into synaptic vesicles before it is released into the synaptic cleft. Mutations in this gene are the cause of autosomal-dominant nonsyndromic type 25 deafness. Alternate splicing results in multiple transcript variants.



## 订购热线: 4008-898-798

Applications Western blotting Predicted band size:65 kDa Positive control:Human cerebella tissue and Human cerebrum tissue lysates Recommended dilution: 200-1000

Gel: 8%SDS-PAGE Lysate: 40 µg Lane 1-2: Human cerebella tissue and Human cerebrum tissue lysates Primary antibody: ml263210(SLC17A8 Antibody) at dilution 1/250 Secondary antibody: Goat anti rabbit IgG at 1/8000 dilution Exposure time: 10 seconds



#### ELISA

Recommended dilution:5000-10000

Wit produce 联系电话: 4008-898-798, 021-61725725

### 联系QQ: 2881505695, 2881505696

- 邮箱: mlbio\_cn@yeah.net
- 网址: www.mlbio.cn