

## Anti-TMEM254 antibody

<b>Cat. No.</b>	ml263313
<b>Package</b>	25 µl/100 µl/200 µl
<b>Storage</b>	-20°C, pH7.4 PBS, 0.05% NaN <sub>3</sub> , 40% Glycerol

### Product overview

<b>Description</b>	Anti-TMEM254 rabbit polyclonal antibody
<b>Applications</b>	ELISA, WB
<b>Immunogen</b>	Synthetic peptide of human TMEM254
<b>Reactivity</b>	Human, Mouse
<b>Content</b>	0.3 mg/ml
<b>Host species</b>	Rabbit
<b>Ig class</b>	Immunogen-specific rabbit IgG
<b>Purification</b>	Antigen affinity purification

### Target information

<b>Symbol</b>	TMEM254
<b>Full name</b>	transmembrane protein 254
<b>Synonyms</b>	C10orf57; bA369J21.6
<b>Swissprot</b>	Q8TBM7

### Target Background

The protein belongs to the TMEM214 family. Subcellular location is Membrane; Multi-pass membrane protein.

订购热线: 4008-898-798

### Applications

#### Western blotting

Predicted band size:14 kDa

Positive control:Human heart tissue, PC-3 cell, Mouse heart tissue and Human cerebrum tissue lysates

Recommended dilution: 200-1000

Gel: 12%SDS-PAGE


Lysate: 40 µg

Lane 1-4: Human heart tissue, PC-3 cell, Mouse heart tissue and Human cerebrum tissue lysates

Primary antibody: ml263313(TM254 Antibody) at dilution 1/300

Secondary antibody: Goat anti rabbit IgG at 1/8000 dilution

Exposure time: 15 seconds



### ELISA

Recommended dilution:5000-10000

联系电话: 4008-898-798, 021-61725725

联系QQ: 2881505695, 2881505696

邮箱: [mlbio\\_cn@yeah.net](mailto:mlbio_cn@yeah.net)

网址: [www.mlbio.cn](http://www.mlbio.cn)