

订购热线: 4008-898-798

Anti-CEP170 antibody

Cat. No.	ml263436
Package	25 μl/100 μl/200 μl
Storage	-20°C, pH7.4 PBS, 0.05% NaN3, 40% Glycerol
Product overview	
Description	Anti-CEP170 rabbit polyclonal antibody
Applications	ELISA, WB, IHC
Immunogen	Synthetic peptide of human CEP170
Reactivity	Human, Mouse
Content	0.5 mg/ml
Host species	Rabbit
lg class	Immunogen-specific rabbit IgG
Purification	Antigen affinity purification
	Antigen affinity purification
Target information	orou
Symbol	CEP170
Full name	centrosomal protein 170
(-,00



Synonyms KAB; FAM68A; KIAA0470

Swissprot Q5SW79

Target Background

The product of this gene is a component of the centrosome, a non-membraneous organelle that functions as the major microtubule-organizing center in animal cells. During interphase, the encoded protein localizes to the sub-distal appendages of mature centrioles, which are microtubule-based structures thought to help organize centrosomes. During mitosis, the protein associates with spindle microtubules near the centrosomes. The protein interacts with and is phosphorylated by polo-like kinase 1, and functions in maintaining microtubule organization and cell morphology. The human genome contains a putative transcribed pseudogene. Several alternatively spliced transcript variants of this gene have been found, but the full-length nature of some of these variants has not been determined.



订购热线: 4008-898-798


Applications

Immunohistochemistry

Predicted cell location: Cytoplasm

Positive control: Human liver cancer

Recommended dilution: 10-50



Good elisa kit producer

The image on the left is immunohistochemistry of paraffin-embedded Human liver cancer tissue using ml263436(CEP170 Antibody) at dilution 1/25, on the right is treated with synthetic peptide. (Original magnification: ×200)

Western blotting

Predicted band size:175 kDa


Positive control:PC-3 cell lysate


Recommended dilution: 200-1000

Predicted cell location: Cytoplasm

Positive control: Human ovarian cancer

Recommended dilution: 10-50





The image on the left is immunohistochemistry of paraffin-embedded Human ovarian cancer tissue using ml263436(CEP170 Antibody) at dilution 1/25, on the right is treated with synthetic peptide. (Original magnification: ×200)



Gel: 6%SDS-PAGE	kDa	
Lysate: 40 µg	250 —	-
Lane: PC-3 cell lysate	130 —	
Primary antibody: ml263436(CEP170 Antibody) at dilution 1/300	95 —	
Secondary antibody: Goat anti rabbit IgG at 1/8000 dilution	72 — 55 — 36 —	
Exposure time: 1 minute		

ELISA

Recommended dilution: 5000-10000

联系电话: 4008-898-798, 021-61725725

联系QQ: 2881505695, 2881505696

邮箱: mlbio_cn@yeah.net

网址: www.mlbio.cn