

订购热线: 4008-898-798

Anti-EPHA6 antibody

Cat. No.	ml263533
Package	25 μl/100 μl/200 μl
Storage	-20°C, pH7.4 PBS, 0.05% NaN3, 40% Glycerol
Product overview	
Description	Anti-EPHA6 rabbit polyclonal antibody
Applications	ELISA, WB, IHC
Immunogen	Synthetic peptide of human EPHA6
Reactivity	Human, Mouse, Rat
Content	0.6 mg/ml
Host species	Rabbit
lg class	Immunogen-specific rabbit IgG
Purification	Antigen affinity purification
	Antigen affinity purification
Target information	oroc
Symbol	EPHA6
Full name	EPH receptor A6
(-,00-



Synonyms EHK2; EK12; EPA6; EHK-2; HEK12; PRO57066

Swissprot Q9UF33

Target Background

Receptor tyrosine kinase which binds promiscuously GPI-anchored ephrin-A family ligands residing on adjacent cells, leading to contact-dependent bidirectional signaling into neighboring cells. The signaling pathway downstream of the receptor is referred to as forward signaling while the signaling pathway downstream of the ephrin ligand is referred to as reverse signaling.



订购热线: 4008-898-798


Applications

Immunohistochemistry

Predicted cell location: Cytoplasm and Cell membrane

Positive control: Human liver cancer

Recommended dilution: 25-100



Good elisakii producer

The image on the left is immunohistochemistry of paraffin-embedded Human liver cancer tissue using ml263533(EPHA6 Antibody) at dilution 1/25, on the right is treated with synthetic peptide. (Original magnification: ×200)

Western blotting

Predicted band size:116 kDa


Positive control:A549 and 293T cell,Mouse brain tissue lysates


Recommended dilution: 500-2000

Predicted cell location: Cytoplasm and Cell membrane

Positive control: Human gastric cancer

Recommended dilution: 25-100





The image on the left is immunohistochemistry of paraffin-embedded Human gastric cancer tissue using ml263533(EPHA6 Antibody) at dilution 1/25, on the right is treated with synthetic peptide. (Original magnification: ×200)



Gel: 6%SDS-PAGE	kDa	1	2	з
Lysate: 40 µg	250 —			
Lane 1-3: A549 and 293T cell, Mouse brain tissue lysates	130—			-
	95 —			
Primary antibody: ml263533(EPHA6 Antibody) at dilution 1/500	72 —			
Secondary antibody: Goat anti rabbit IgG at 1/8000 dilution	55 —			
	36 —	—		
Exposure time: 10 seconds				

ELISA

Recommended dilution: 5000-10000

联系电话: 4008-898-798, 021-61725725

联系QQ: 2881505695, 2881505696

邮箱: mlbio_cn@yeah.net

网址: www.mlbio.cn