

Anti-SLC5A8 antibody

 Cat. No.
 ml263602

 Package
 25 μl/100 μl/200 μl

 Storage
 -20°C, pH7.4 PBS, 0.05% NaN3, 40% Glycerol

Product overview	
Description	Anti-SLC5A8 rab
Applications	ELISA, WB
Immunogen	Synthetic peptide
Reactivity	Human
Content	0.7 mg/ml
Host species	Rabbit
lg class	Immunogen-spec
Purification	Antigen affinity p
Target information	

-100de

Anti-SLC5A8 rabbit polyclonal antibody ELISA, WB Synthetic peptide of human SLC5A8 Human D.7 mg/ml Rabbit mmunogen-specific rabbit IgG Antigen affinity purification

Target information Symbol Full name Synonyms Swissprot SLC5A8 solute carrier family 5 member 8 AIT; SMCT; SMCT1 Q8N695

Target Background

SLC5A8 has been shown to transport iodide by a passive mechanism and monocarboxylates and short-chain fatty acids by a sodium-coupled mechanism .In kidney, SLC5A8 functions as a high-affinity sodium-coupled lactate transporter involved in reabsorption of lactate and maintenance of blood lactate levels.



订购热线: 4008-898-798

Applications Western blotting Predicted band size:67 kDa Positive control:Human thyroid tissue lysate Recommended dilution: 500-2000

Gel: 8%SDS-PAGE Lysate: 40 µg Lane: Human thyroid tissue lysate Primary antibody: ml263602(SLC5A8 Antibody) at dilution 1/600 Secondary antibody: Goat anti rabbit IgG at 1/8000 dilution Exposure time: 40 seconds

28

ELISA

Recommended dilution:5000-10000

Wit produce 联系电话: 4008-898-798, 021-61725725

联系QQ: 2881505695, 2881505696

- 邮箱: mlbio_cn@yeah.net
- 网址: www.mlbio.cn