

Anti-ACD antibody

 Cat. No.
 ml263606

 Package
 25 μl/100 μl/200 μl

 Storage
 -20°C, pH7.4 PBS, 0.05% NaN3, 40% Glycerol

Product overview	
Description	Anti-ACD rabbit polyclonal antibody
Applications	ELISA, WB
Immunogen	Synthetic peptide of human ACD
Reactivity	Human, Mouse, Rat
Content	1 mg/ml
Host species	Rabbit
lg class	Immunogen-specific rabbit IgG
Purification	Antigen affinity purification
Target information	
Symbol	ACD
Full name	adrenocortical dysplasia homolog

adrenocortical dysplasia homolo PIP1; PTOP; TPP1; TINT1 Q96AP0

Target Background

Synonyms

Swissprot


This gene encodes a protein that is involved in telomere function. This protein is one of six core proteins in the telosome/shelterin telomeric complex, which functions to maintain telomere length and to protect telomere ends. Through its interaction with other components, this protein plays a key role in the assembly and stabilization of this complex, and it mediates the access of telomerase to the telomere. Multiple transcript variants encoding different isoforms have been found for this gene. This gene, which is also referred to as TPP1, is distinct from the unrelated TPP1 gene on chromosome 11, which encodes tripeptidyl-peptidase I.



订购热线: 4008-898-798

Applications Western blotting Predicted band size:58 kDa Positive control:HEPG2 cell lysate Recommended dilution: 500-2000

Gel: 8%SDS-PAGE Lysate: 40 µg Lane: HEPG2 cell lysate Primary antibody: ml263606(ACD Antibody) at dilution 1/400 Secondary antibody: Goat anti rabbit IgG at 1/8000 dilution Exposure time: 60 seconds



ELISA

Recommended dilution:5000-10000

Wit produces 联系电话: 4008-898-798, 021-61725725

联系QQ: 2881505695, 2881505696

- 邮箱: mlbio_cn@yeah.net
- 网址: www.mlbio.cn