

Anti-CDX4 antibody

 Cat. No.
 ml263618

 Package
 25 μl/100 μl/200 μl

 Storage
 -20°C, pH7.4 PBS, 0.05% NaN3, 40% Glycerol

Product overview	
Description	Anti-CDX4 rabbit polyclonal antibody
Applications	ELISA, WB
Immunogen	Synthetic peptide of human CDX4
Reactivity	Human, Mouse
Content	0.3 mg/ml
Host species	Rabbit
lg class	Immunogen-specific rabbit IgG
Purification	Antigen affinity purification
Target information	
Symbol	CDX4
Full name	caudal type homeobox 4
Synonyms	

O14627

Target Background

Swissprot

Cdx1, Cdx2 and Cdx4 are members of the caudal-type homeobox family of genes, which are homologs of the Drosophila "caudal" gene required for anterior-posterior regional identity. The proteins encoded by these genes are transcription factors which play an important role in development by regulating the expression of Hox genes. Hox genes play a fundamental role in the development of the vertebrate central nervous system, heart, axial skeleton, limbs, gut, urogenital tract and external genitalia. Cdx4 is a major positive regulator of the expression of all Hox family members. Due to its critical role as a regulator, Cdx4 is a direct target of the canonical Wnt pathway. The loss of Cdx4 can result in the development of an expanded hindbrain, while the overexpression of Cdx4 may cause the hindbrain to lose its distinct segmental features and resemble the spinal cord.



订购热线: 4008-898-798

Applications Western blotting Predicted band size:30 kDa Positive control:Hela and HEPG2 cell, Human fetal brain tissue lysates Recommended dilution: 200-1000

kDa 1 2 з 130 Gel: 8%SDS-PAGE Lysate: 40 µg 95 -Lane 1-3: Hela and HEPG2 cell, Human fetal brain tissue 72 lysates 55 Primary antibody: ml263618(CDX4 Antibody) at dilution 1/200 36 Secondary antibody: Goat anti rabbit IgG at 1/8000 dilution 28 Exposure time: 5 minutes 17

ELISA

Recommended dilution:5000-10000

联系电话: 4008-898-798, 021-61725725

联系QQ: 2881505695, 2881505696

- 邮箱: mlbio_cn@yeah.net
- 网址: www.mlbio.cn