

订购热线: 4008-898-798

# Anti-COX6C antibody

**Cat. No.** ml263662

**Package** 25 μl/100 μl/200 μl

**Storage** -20°C, pH7.4 PBS, 0.05% NaN3, 40% Glycerol

**Product overview** 

**Description** Anti-COX6C rabbit polyclonal antibody

Applications ELISA, WB

Immunogen Synthetic peptide of human COX6C

ReactivityHuman, MouseContent0.7 mg/mlHost speciesRabbit

Ig classImmunogen-specific rabbit IgGPurificationAntigen affinity purification

**Target information** 

Symbol COX6C

Full name cytochrome c oxidase subunit 6C

Synonyms

Swissprot P09669

## **Target Background**

Cytochrome c oxidase, the terminal enzyme of the mitochondrial respiratory chain, catalyzes the electron transfer from reduced cytochrome c to oxygen. It is a heteromeric complex consisting of 3 catalytic subunits encoded by mitochondrial genes and multiple structural subunits encoded by nuclear genes. The mitochondrially-encoded subunits function in electron transfer, and the nuclear-encoded subunits may be involved in the regulation and assembly of the complex. This nuclear gene encodes subunit VIc, which has 77% amino acid sequence identity with mouse subunit VIc. This gene is upregulated in prostate cancer cells. A pseudogene has been found on chromosomes 16p12.



订购热线: 4008-898-798

## **Applications**

## Western blotting

Predicted band size:9 kDa

Positive control:Mouse heart tissue and Human heart tissue lysates

Recommended dilution: 500-2000

Gel: 12%SDS-PAGE

Lysate: 40 µg

Lane 1-2: Mouse heart tissue and Human heart tissue


lysates

Primary antibody: ml263662(COX6C Antibody) at dilution

1/800

Secondary antibody: Goat anti rabbit IgG at 1/8000 dilution

Exposure time: 1 second



## **ELISA**

Recommended dilution:5000-10000

联系电话: 4008-898-798, 021-61725725

联系QQ: 2881505695, 2881505696

邮箱: mlbio\_cn@yeah.net

网址: www.mlbio.cn