

Anti-DZIP1 antibody

Cat. No. ml263751

Package 25 μl/100 μl/200 μl

Storage -20°C, pH7.4 PBS, 0.05% NaN3, 40% Glycerol

Product overview

Description Anti-DZIP1 rabbit polyclonal antibody

Applications ELISA, WB, IHC

Immunogen Synthetic peptide of human DZIP1

Reactivity Human

Content 1.56 mg/ml

Host species Rabbit

Ig class Immunogen-specific rabbit IgG

Purification Antigen affinity purification

Target information

Symbol DZIP1

Full name DAZ interacting zinc finger protein 1



Synonyms DZIP; DZIPt1

Swissprot Q86YF9

Target Background

DZIP1(DAZ interacting protein 1) is also known as DZIP or DZIP2 and is a 867 amino acid protein which is expressed as 3 isoforms, designated DZIPb, DZIPt1 and DZIPt2. DZIP1 is localized to testis, oocytes, ovary and fetal brain, as well as in embryonic stem cells and germ cells. In testis, DZIP1 is localized to the nucleus and also shows some cytoplasmic distribution in spermatogonia. DZIP1 belongs to the C2H2-type zinc-finger protein family, and, characteristic of the C2H2-type zinc-finger protein family, DZIP1 contains one C2H2-type zinc-finger region through which it is thought to interact with DAZ, an interaction that promotes spermiogenesis. DZIP1 expression is not found in those afflicted with Sertoli cell-only syndrome (characterized by the absence of germ cells in the testis), suggesting that the lack of DZIP1 may be involved in the pathogenesis of Sertoli cell-only syndrome-induced male sterility




Applications


Immunohistochemistry

Predicted cell location: Cytoplasm


Positive control: Human liver cancer

Recommended dilution: 25-100





Good elisakii produceri



Western blotting

Predicted band size:99 kDa

Positive control:PC-3 cell lysate

Recommended dilution: 200-1000

Predicted cell location: Cytoplasm

Positive control: Human brain

Recommended dilution: 25-100





The image on the left is immunohistochemistry of paraffin-embedded Human brain tissue using ml263751(DZIP1 Antibody) at dilution 1/40, on the right is treated with synthetic peptide. (Original magnification: ×200)



Gel: 6%SDS-PAGE


Lysate: 40 µg

Lane: PC-3 cell lysate

Primary antibody: ml263751(DZIP1 Antibody) at dilution 1/200

Secondary antibody: Goat anti rabbit IgG at 1/8000 dilution

Exposure time: 5 minutes



ELISA

Recommended dilution: 5000-10000

联系电话: 4008-898-798, 021-61725725

联系QQ: 2881505695, 2881505696

邮箱: mlbio_cn@yeah.net

网址: www.mlbio.cn