

Anti-MYDGF antibody

Cat. No.	ml263807
Package	25 µl/100 µl/200 µl
Storage	-20°C, pH7.4 PBS, 0.05% NaN ₃ , 40% Glycerol

Product overview

Description	Anti-MYDGF rabbit polyclonal antibody
Applications	ELISA, WB
Immunogen	Synthetic peptide of human MYDGF
Reactivity	Human, Mouse
Content	0.66 mg/ml
Host species	Rabbit
Ig class	Immunogen-specific rabbit IgG
Purification	Antigen affinity purification

Target information

Symbol	MYDGF
Full name	myeloid derived growth factor
Synonyms	IL25; IL27; SF20; IL27w; C19orf10; R33729_1; EUROIMAGE1875335
Swissprot	Q969H8

Target Background

The protein encoded by this gene was previously thought to support proliferation of lymphoid cells and was considered an interleukin. However, this activity has not been reproducible and the function of this protein is currently unknown.

订购热线: 4008-898-798

Applications

Western blotting

Predicted band size:19 kDa

Positive control:HEPG2 cell, Human between peritoneal stromal sarcoma tissue, K562 and 293T cell lysates

Recommended dilution: 500-2000

Gel: 12%SDS-PAGE


Lysate: 40 µg

Lane 1-4: HEPG2 cell, Human between peritoneal stromal sarcoma tissue, K562 and 293T cell lysates

Primary antibody: ml263807(MYDGF Antibody) at dilution 1/250

Secondary antibody: Goat anti rabbit IgG at 1/8000 dilution

Exposure time: 6 minutes



ELISA

Recommended dilution:5000-10000

联系电话: 4008-898-798, 021-61725725

联系QQ: 2881505695, 2881505696

邮箱: mlbio_cn@yeah.net

网址: www.mlbio.cn

mlbio 酶联生物
Good elisakit producers
订购热线: 4008-898-798