

Anti-FKTN antibody

Cat. No.	ml263948
Package	25 μl/100 μl/200 μl
Storage	-20°C, pH7.4 PBS, 0.05% NaN3, 40% Glycerol

Product overview	
Description	Anti-FKTN rabbit polyclonal antibody
Applications	ELISA, WB
Immunogen	Synthetic peptide of human FKTN
Reactivity	Human, Mouse
Content	1.08 mg/ml
Host species	Rabbit
lg class	Immunogen-specific rabbit IgG
Purification	Antigen affinity purification
Target information	
Symbol	FKTN
Full name	fukutin
Synonyms	FCMD; CMD1X; LGMD2M; MDDGA4; MDDGB4; MDDGC4
Swissprot	O75072
	C ^C
Target Background	

Target Background

The protein encoded by this gene is a putative transmembrane protein that is localized to the cis-Golgi compartment, where it may be involved in the glycosylation of alpha-dystroglycan in skeletal muscle. The encoded protein is thought to be a glycosyltransferase and could play a role in brain development. Defects in this gene are a cause of Fukuyama-type congenital muscular dystrophy (FCMD), Walker-Warburg syndrome (WWS), limb-girdle muscular dystrophy type 2M (LGMD2M), and dilated cardiomyopathy type 1X (CMD1X). Alternatively spliced transcript variants have been found for this gene.



订购热线: 4008-898-798

Applications Western blotting Predicted band size:54 kDa Positive control:Mouse liver tissue lysate Recommended dilution: 500-2000

Gel: 8%SDS-PAGE Lysate: 40 µg Lane: Mouse liver tissue lysate Primary antibody: ml263948(FKTN Antibody) at dilution 1/450 Secondary antibody: Goat anti rabbit IgG at 1/8000 dilution Exposure time: 1 minute

28

17



ELISA

Recommended dilution:5000-10000

Wit produces 联系电话: 4008-898-798, 021-61725725

联系QQ: 2881505695, 2881505696

- 邮箱: mlbio_cn@yeah.net
- 网址: www.mlbio.cn