

Anti-ZNF85 antibody

Cat. No. ml264157
Package 25 µl/100 µl/200 µl
Storage -20°C, pH7.4 PBS, 0.05% NaN₃, 40% Glycerol

Product overview

Description Anti-ZNF85 rabbit polyclonal antibody
Applications ELISA, WB
Immunogen Synthetic peptide of human ZNF85
Reactivity Human, Mouse
Content 0.84 mg/ml
Host species Rabbit
Ig class Immunogen-specific rabbit IgG
Purification Antigen affinity purification

Target information

Symbol ZNF85
Full name zinc finger protein 85
Synonyms HPF4; HTF1
Swissprot Q03923

Target Background

May be a transcriptional repressor.

订购热线: 4008-898-798

Applications

Western blotting

Predicted band size:69 kDa

Positive control:Human fetal liver tissue, PC3, Raji, HeLa and TM4 cell lysates

Recommended dilution: 500-2000

Gel: 6%SDS-PAGE


Lysate: 40 µg

Lane 1-5: Human fetal liver tissue, PC3, Raji, HeLa and TM4 cell lysates

Primary antibody: ml264157(ZNF85 Antibody) at dilution 1/350

Secondary antibody: Goat anti rabbit IgG at 1/8000 dilution

Exposure time: 10 seconds



ELISA

Recommended dilution:5000-10000

联系电话: 4008-898-798, 021-61725725

联系QQ: 2881505695, 2881505696

邮箱: mlbio_cn@yeah.net

网址: www.mlbio.cn