

订购热线: 4008-898-798

Anti-HES7 antibody

Cat. No.	ml264171
Package	25 μl/100 μl/200 μl
Storage	-20°C, pH7.4 PBS, 0.05% NaN3, 40% Glycerol
Product overview	
Description	Anti-HES7 rabbit polyclonal antibody
Applications	ELISA, WB, IHC
Immunogen	Synthetic peptide of human HES7
Reactivity	Human, Mouse
Content	0.72 mg/ml
Host species	Rabbit
lg class	Immunogen-specific rabbit IgG
Purification	Antigen affinity purification
	Antigen affinity purification
Target information	oroc
Symbol	HES7
Full name	hes family bHLH transcription factor 7
C	,0 ⁰ ¹



Synonyms SCDO4; bHLHb37

Swissprot Q9BYE0

Target Background

This gene encodes a member of the hairy and enhancer of split family of bHLH transcription factors. The mouse ortholog of this gene is regulated by Notch signaling. The protein functions as a transcriptional repressor, and is implicated in correct patterning of the axial skeleton. A mutation in this gene has been shown to result in spondylocostal dysostosis. Multiple transcript variants encoding different isoforms have been found for this gene.



订购热线: 4008-898-798


Applications

Immunohistochemistry

Predicted cell location: Nucleus

Positive control: Human liver cancer

Recommended dilution: 30-150



The image on the left is immunohistochemistry of paraffin-embedded Human liver cancer tissue using ml264171(HES7 Antibody) at dilution 1/25, on the right is treated with synthetic peptide. (Original magnification: ×200)

Western blotting

Predicted band size:25 kDa

Positive control:Hela cell lysate

Recommended dilution: 200-1000

Good elisa Nil producers



Gel: 12%SDS-PAGE	kDa	
Lysate: 40 µg	72— 55—	
	36—	
Lane: Hela cell lysate	28—	
Primary antibody: ml264171(HES7 Antibody) at dilution 1/250	17—	-
Secondary antibody: Goat anti rabbit IgG at 1/8000 dilution		
Exposure time: 1 minute	10—	

ELISA

Recommended dilution: 5000-10000

联系电话: 4008-898-798, 021-61725725

联系QQ: 2881505695, 2881505696

邮箱: mlbio_cn@yeah.net

网址: www.mlbio.cn