

Anti-HNRNPA1L2 antibody

Cat. No. ml264172

Package 25 μl/100 μl/200 μl

Storage -20°C, pH7.4 PBS, 0.05% NaN3, 40% Glycerol

Product overview

Description Anti-HNRNPA1L2 rabbit polyclonal antibody

Applications ELISA, WB, IHC

Immunogen Synthetic peptide of human HNRNPA1L2

Reactivity Human

Content 0.6 mg/ml

Host species Rabbit

Ig class Immunogen-specific rabbit IgG

Purification Antigen affinity purification

Target information

Symbol HNRNPA1L2

Full name heterogeneous nuclear ribonucleoprotein A1-like 2



S	ν	n	o	n	ν	m	าร

Swissprot Q32P51

Target Background

Involved in the packaging of pre-mRNA into hnRNP particles, transport of poly(A) mRNA from the nucleus to the cytoplasm and may modulate splice site selection.




Applications


Immunohistochemistry


Predicted cell location: Nucleus Predicted cell location: Nucleus


Positive control: Human tonsil Positive control: Human breast cancer

Recommended dilution: 20-100 Recommended dilution: 20-100









The image on the left is immunohistochemistry of paraffin-embedded Human tonsil tissue using ml264172(HNRNPA1L2 Antibody) at dilution 1/20, on the right is treated with synthetic peptide. (Original magnification: ×200)

The image on the left is immunohistochemistry of paraffin-embedded Human breast cancer tissue using ml264172(HNRNPA1L2 Antibody) at dilution 1/20, on the right is treated with synthetic peptide. (Original magnification: ×200)

Western blotting

Predicted band size:34 kDa

Positive control:293T, K562, HepG2, 231, Hela and Jurkat cell lysates

Good elisakit producer

Recommended dilution: 500-2000



Gel: 8%SDS-PAGE

Lysate: 40 µg


Lane 1-6: 293T, K562, HepG2, 231, Hela and Jurkat cell lysates

Primary antibody: ml264172(HNRNPA1L2 Antibody) at dilution

1/250

Secondary antibody: Goat anti rabbit IgG at 1/8000 dilution

Exposure time: 1 minute



ELISA

Recommended dilution: 5000-10000

联系电话: 4008-898-798, 021-61725725

联系QQ: 2881505695, 2881505696

邮箱: mlbio_cn@yeah.net

网址: www.mlbio.cn