

Anti-HPDL antibody

Cat. No.	ml264214
Package	25 µl/100 µl/200 µl
Storage	-20°C, pH7.4 PBS, 0.05% NaN ₃ , 40% Glycerol

Product overview

Description	Anti-HPDL rabbit polyclonal antibody
Applications	ELISA, WB
Immunogen	Synthetic peptide of human HPDL
Reactivity	Human
Content	0.9 mg/ml
Host species	Rabbit
Ig class	Immunogen-specific rabbit IgG
Purification	Antigen affinity purification

Target information

Symbol	HPDL
Full name	4-hydroxyphenylpyruvate dioxygenase like
Synonyms	GLOXD1; 4-HPPD-L
Swissprot	Q96IR7

Target Background

May have dioxygenase activity.

订购热线: 4008-898-798

Applications

Western blotting

Predicted band size:39 kDa

Positive control:Hela and 231 cell lysates

Recommended dilution: 500-2000

Gel: 8%SDS-PAGE


Lysate: 40 μ g

Lane 1-2: Hela and 231 cell lysates

Primary antibody: ml264214(HPDL Antibody) at dilution 1/400

Secondary antibody: Goat anti rabbit IgG at 1/8000 dilution

Exposure time: 1 minute



ELISA

Recommended dilution:5000-10000

联系电话: 4008-898-798, 021-61725725

联系QQ: 2881505695, 2881505696

邮箱: mlbio_cn@yeah.net

网址: www.mlbio.cn