

订购热线: 4008-898-798

# Anti-RRP8 antibody

**Cat. No.** ml264224

Package 25  $\mu$ l/100  $\mu$ l/200  $\mu$ l

**Storage** -20°C, pH7.4 PBS, 0.05% NaN3, 40% Glycerol

**Product overview** 

**Description** Anti-RRP8 rabbit polyclonal antibody

Applications ELISA, WB

**Immunogen** Synthetic peptide of human RRP8

Reactivity Human, Mouse, Rat

Content 1.02 mg/ml Host species Rabbit

Ig classImmunogen-specific rabbit IgGPurificationAntigen affinity purification

**Target information** 

Symbol RRP8

Full name ribosomal RNA processing 8

Synonyms NML; KIAA0409

Swissprot O43159

# **Target Background**

Essential component of the eNoSC (energy-dependent nucleolar silencing) complex, a complex that mediates silencing of rDNA in response to intracellular energy status and acts by recruiting histone-modifying enzymes. The eNoSC complex is able to sense the energy status of cell: upon glucose starvation, elevation of NAD+/NADP+ ratio activates SIRT1, leading to histone H3 deacetylation followed by dimethylation of H3 at 'Lys-9' (H3K9me2) by SUV39H1 and the formation of silent chromatin in the rDNA locus. In the complex, RRP8 binds to H3K9me2 and probably acts as a methyltransferase. Its substrates are however unknown.



订购热线: 4008-898-798

# **Applications**

# Western blotting

Predicted band size:51 kDa

Positive control:Human bladder transitional cell carcinoma grade 2-3 tissue lysate

Recommended dilution: 1000-5000

Gel: 8%SDS-PAGE

Lysate: 40 µg

Lane: Human bladder transitional cell carcinoma grade 2-3

tissue lysate

Primary antibody: ml264224(RRP8 Antibody) at dilution

1/800

Secondary antibody: Goat anti rabbit IgG at 1/5000 dilution

Exposure time: 1 minute

kDa 250— 130— 95— 72— 55— 36—

# **ELISA**

Recommended dilution:5000-10000

联系电话: 4008-898-798, 021-61725725

联系QQ: 2881505695, 2881505696

邮箱: mlbio\_cn@yeah.net

网址: www.mlbio.cn