

订购热线: 4008-898-798

Anti-IRX1 antibody

Cat. No.	ml264274
Package	25 µl/100 µl/200 µl
Storage	-20°C, pH7.4 PBS, 0.05% NaN3, 40% Glycerol
Product overview	
Description	Anti-IRX1 rabbit polyclonal antibody
Applications	ELISA, WB, IHC
Immunogen	Synthetic peptide of human IRX1
Reactivity	Human, Mouse
Content	1.2 mg/ml
Host species	Rabbit
lg class	Immunogen-specific rabbit IgG
Purification	Antigen affinity purification
	Antigen attinity purification
Target information	orot
Symbol	IRX1
Full name	iroquois homeobox 1
(,00 [°]



Synonyms IRX-5; IRXA1

Swissprot P78414

Target Background

This gene encodes a member of the Iroquois homeobox protein family. Homeobox genes in this family are involved in pattern formation in the embryo. The gene product has been identified as a tumor suppressor in gastric and head and neck cancers. A pseudogene of this gene is located on chromosome 13.



订购热线: 4008-898-798


Applications

Immunohistochemistry

Predicted cell location: Nucleus

Positive control: Human esophagus cancer

Recommended dilution: 40-200



Good elisa Nit producert

The image on the left is immunohistochemistry of paraffin-embedded Human esophagus cancer tissue using ml264274(IRX1 Antibody) at dilution 1/40, on the right is treated with synthetic peptide. (Original magnification: ×200)

Western blotting

Predicted band size:50 kDa


Positive control:Human heart tissue lysate


Recommended dilution: 500-2000

Predicted cell location: Nucleus

Positive control: Human ovarian cancer

Recommended dilution: 40-200





The image on the left is immunohistochemistry of paraffin-embedded Human ovarian cancer tissue using ml264274(IRX1 Antibody) at dilution 1/40, on the right is treated with synthetic peptide. (Original magnification: ×200)



Gel: 8%SDS-PAGE

ELISA

Recommended dilution: 5000-10000

联系电话: 4008-898-798,021-61725725

联系QQ: 2881505695, 2881505696

邮箱: mlbio_cn@yeah.net

网址: www.mlbio.cn