

订购热线: 4008-898-798

# Anti-PPY antibody

**Cat. No.** ml264297

**Package** 25 μl/100 μl/200 μl

**Storage** -20°C, pH7.4 PBS, 0.05% NaN3, 40% Glycerol

**Product overview** 

**Description** Anti-PPY rabbit polyclonal antibody

Applications ELISA, WB

**Immunogen** Synthetic peptide of human PPY

ReactivityHuman, RatContent2.28 mg/mlHost speciesRabbit

Ig classImmunogen-specific rabbit IgGPurificationAntigen affinity purification

**Target information** 

Symbol PPY

Full name pancreatic polypeptide

SynonymsPP; PNPSwissprotP01298

### **Target Background**

This gene encodes a member of the neuropeptide Y (NPY) family of peptides. The encoded 95 aa preproprotein is synthesized in the pancreatic islets of Langerhans and proteolytically processed to generate two peptide products. These products include the active pancreatic hormone of 36 aa and an icosapeptide of unknown function. This hormone acts as a regulator of pancreatic and gastrointestinal functions and may be important in the regulation of food intake. Plasma level of this hormone has been shown to be reduced in conditions associated with increased food intake and elevated in anorexia nervosa. In addition, infusion of this hormone in obese rodents has shown to decrease weight gain. Alternative splicing results in multiple transcript variants, at least one of which encodes an isoform that is proteolytically processed.



订购热线: 4008-898-798

## **Applications**

## Western blotting

Predicted band size:11 kDa

Positive control:Rat pancreas tissue lysate

Recommended dilution: 1000-5000

Gel: 12%SDS-PAGE

Lysate: 40 µg

Lane: Rat pancreas tissue lysate

Primary antibody: ml264297(PPY Antibody) at dilution

1/1900

Secondary antibody: Goat anti rabbit IgG at 1/5000 dilution

Exposure time: 3 minutes

kDa
95—
72—
55—
36—
28—
17—
10—

### **ELISA**

Recommended dilution:5000-10000

联系电话: 4008-898-798, 021-61725725

联系QQ: 2881505695, 2881505696

邮箱: mlbio\_cn@yeah.net

网址: www.mlbio.cn