

订购热线: 4008-898-798

Anti-JRKL antibody

Cat. No. ml264302

Package 25 μl/100 μl/200 μl

Storage -20°C, pH7.4 PBS, 0.05% NaN3, 40% Glycerol

Product overview

Description Anti-JRKL rabbit polyclonal antibody

Applications ELISA, WB

Immunogen Synthetic peptide of human JRKL

ReactivityHuman, MouseContent1.38 mg/mlHost speciesRabbit

Ig classImmunogen-specific rabbit IgGPurificationAntigen affinity purification

Target information

Symbol JRKL
Full name JRK like
Synonyms HHMJG
Swissprot Q9Y4A0

Target Background

The function of this gene has not yet been defined, however, the encoded protein shares similarity with the human (41% identical) and mouse (34% identical) jerky gene products. This protein may act as a nuclear regulatory protein. Alternatively spliced transcript variants have been found for this gene.



订购热线: 4008-898-798

Applications

Western blotting

Predicted band size:60 kDa
Positive control:K562 cell lysate
Recommended dilution: 500-2000

Gel: 8%SDS-PAGE

Lysate: 40 µg

Lane: K562 cell lysate

Primary antibody: ml264302(JRKL Antibody) at dilution

1/1000

Secondary antibody: Goat anti rabbit IgG at 1/5000 dilution

Exposure time: 7 minutes

kDa 250— 130— 95— 72— 55—

ELISA

Recommended dilution:5000-10000

联系电话: 4008-898-798, 021-61725725

联系QQ: 2881505695, 2881505696

邮箱: mlbio_cn@yeah.net

网址: www.mlbio.cn