

高范围 RNA Ladder 说明书

中文名称：高范围 RNA Ladder

英文名称：RiboRuler High range RNA Ladder

级别：BR

分子量范围：200~6000

产品描述： The RiboRuler High range RNA Ladder is a mixture of eight chromatography-purified single-stranded RNA transcripts. The RNA concentration of each ladder band is precisely measured by spectrophotometry, which makes this ladder ideal for approximate RNA

quantification on gels. The RiboRuler™ High range RNA Ladder is recommended for electrophoresis in the following: native 1% agarose with TAE buffer, denaturing formaldehyde agarose with MOPS buffer, denaturing glyoxal/DMSO agarose with sodium phosphate buffer

The ladder is supplied with 2× RNA Loading Dye solution

Fragments(in bases): 6000, 4000, 3000, 2000, 1500, 1000, 500, 200

性状(以下信息仅供参考): 液体

用途: 本品仅供科研, 不得用于其它用途

保存: -20° C