

氨基琼脂糖凝胶 4B 说明书

中文名称： 氨基琼脂糖凝胶 4B

英文名称： EAH Sepharose 4B

级别： BR

Active group: amino group

Active group concentration: 7-12umoles amino groups/ml drained matrix

Spacer: 1.6-diaminohexane(11 atom)

Coupling capacity: Application dependent

Bead structure: 4% agarose

Bead size range: 45~165um

Mean bead size: 90um

Max linear flow rate: 75cm/h at 25°C; HR 16/10 column; 5cm bed height

pH stability: long term and working, short term : 3~14

Chemical stability: Stable to all commonly used aqueous buffers

Physical stability: Negligible volume variation due to changes in pH or ionic strength



性状(以下信息仅供参考): 含 20%乙醇，底部为乳白色胶体

用途: 本品仅供科研，不得用于其它用途。(以下用途仅供参考)

保存: 2~8℃