

丙烯葡聚糖凝胶 S-1000 SF 说明书

中文名称：丙烯葡聚糖凝胶 S-1000 SF

英文名称：Sephacryl S-1000 Superfine

级别：BR

Principle: Gel filtration

Matrix: Spherical allyl dextran and N,N-methylenebisacrylamide

Column buffer: Distilled water containing 20% ethanol as a preservative

Fractionation range for dextrans(Mr): 500000~100000000



Exclusion limit for DNA: 20000 base pairs

Particle size range: 40~105um

Maximum operating flow rate(approx): 40cm/h

性状(以下信息仅供参考): 含 20%乙醇，底部为胶体。

用途: 本品仅供科研，不得用于其它用途。(以下用途仅供参考)

保存: 2~8℃