

3-磷酸甘油酸激酶说明书

中文名称：3-磷酸甘油酸激酶

英文名称：3-Phosphoglyceric Phosphokinase from baker's yeast ; 3-Phosphoglycerate kinase; ATP:3-Phospho-D-glycerate 1-phosphotransferase

CAS 号：9001-83-6

级别：BR

活力： ≥ 1000 units/mg protein

活力定义：One unit will convert 1.0μ mole of 1,3-diphosphoglycerate

to 3-phosphoglycerate per min at pH 6.9 at 25 °C

性状(以下信息仅供参考): 白色或棕褐色悬浮液, Crystalline suspension in 3.0 M (NH₄)₂SO₄ and 0.04 M tetrasodium pyrophosphate solution, pH 8.0

用途: 本品仅供科研, 不得用于其它用途

保存: 2~8 °C