

甘油醛-3-磷酸脱氢酶说明书

中文名称：甘油醛-3-磷酸脱氢酶

英文名称：GAPDH; Glyceraldehyde-3-phosphate Dehydrogenase from rabbit muscle; D-Glyceraldehyde 3-phosphate:NAD⁺ oxidoreductase (phosphorylating)

其他名称：3-磷酸甘油醛脱氢酶

CAS 号：9001-50-7

级别：BR

活力：≥75units/mg protein

活力定义: One unit will reduce 1.0 μ mole of 3-phosphoglycerate to D-glyceraldehyde 3-phosphate per min in a coupled system with 3-phosphoglyceric phosphokinase at pH 7.6 at 25 °C

Protein: \geq 80% biuret

性状(以下信息仅供参考): 白色或类白色粉末

用途: 本品仅供科研, 不得用于其它用途。(以下用途仅供参考)

保存: -20°C