

HindIII 说明书

中文名称: HindIII

英文名称: HindIII

其他名称: 限制性内切酶 HindIII

级别: BR

来源: Recombinant E. coli strain that carries the cloned hindIII gene from Haemophilus influenzae.

保存条件: 10 mM Tris-HCl (pH 7.5 at 25° C), 250 mM KCl, 1 mM DTT, 0.1 mM EDTA, 0.2 mg/mL BSA and 50% glycerol.

缓冲液 Optimal buffer(10X): 100 mM Tris-HCl (pH 8.5), 100 mM MgCl₂, 1000 mM KCl, 1 mg/mL BSA..

Ligations: After 10-fold over digestion with enzyme 90% of the DNA fragments can be ligated and recut.

使用说明: To obtain 100% activity, BSA should be added to the 1 x reaction mix to a final concentration of 100 μg/ml. Do not use BSA for long incubation.

性状(以下信息仅供参考): 液体

用途: 本品仅供科研，不得用于其它用途。

保存: -20℃