

# 宝石红蛋白染色试剂说明书

中文名称：宝石红蛋白染色试剂

英文名称：SYPRO red protein stain

级别：BS

产品特点：

**1.Highly sensitive** : The stain can detect 1 to 2 ng of protein per mini gel band, making it more sensitive than Coomassie Brilliant Blue or silver staining

**2.Rapid:** Staining is complete in less than one hour

**3.Simple:** Selective for proteins

4. SYPRO Red detects proteins as small as 6.5kDa and does not stain nucleic acids or lipopolysaccharides. Broad linear range of detection

性状(以下信息仅供参考)液体

**用途:** 本品仅供科研，不得用于其它用途。

**保存:** 2~8℃