

双向电泳指示剂说明书

中文名称：双向电泳指示剂

英文名称：Tracking Dye-2D

级别：Proteomics

折光率：1.37~1.41

比重：1.149~1.159

Protease: None detected



产品描述: Premixed,ready-to-use,and a great complimentary product.This tracking dye very accurately monitors the progress of the 2nd dimension gel during the 2D electrophoresis process.

性状(以下信息仅供参考): 澄清紫色液体

用途: 本品仅供科研，不得用于其它用途。

保存: -20℃