

## 3-磷酸甘油醛脱氢酶抗体

产品货号: mIR188

英文名称: GAPDH-Loading Control

中文名称: 3-磷酸甘油醛脱氢酶抗体

别 名: 38 kDa BFA-dependent ADP-ribosylation substrate; Aging-associated gene 9 protein; BARS-38; cb609; EC 1.2.1.12; G3PD; G3PDH; GAPD; Glyceraldehyde 3 phosphate dehydrogenase; Glyceraldehyde 3 phosphate dehydrogenase liver; Glyceraldehyde 3 phosphate dehydrogenase muscle; KNC-NDS6; MGC102544; MGC102546; MGC103190; MGC103191; MGC105239; MGC127711; MGC88685; OCAS, p38 component; OCT1 coactivator in S phase, 38-KD component; wu:fb33a10.

产品类型: 内参抗体

研究领域: 免疫学

抗体来源: Rabbit

克隆类型: Polyclonal

交叉反应 : Human, Mouse, Rat, Rabbit,



**产品应用:** WB=1:2000-10000 IHC-P=1:400-800 IHC-F=1:400-800 IF=1:100-500 (石蜡切片需做抗原修复) not yet tested in other applications.

optimal dilutions/concentrations should be determined by the end user.

分子量: 37kDa

细胞定位: 细胞核 细胞浆 细胞膜

性 状: Lyophilized or Liquid

浓 度: 1mg/ml

免疫原: GAPDH protein of rabbit:

亚 型: IgG

纯化方法: affinity purified by Protein A

储存液: 0.01M TBS(pH7.4) with 1% BSA, 0.03% Proclin300 and 50% Glycerol.

保存条件: Store at -20  $^{\circ}$  C for one year. Avoid repeated freeze/thaw cycles. The lyophilized antibody is stable at room temperature for at least one month and for greater than a year when kept at -20 $^{\circ}$  C. When reconstituted in sterile pH 7.4 0.01M PBS or diluent of antibody the antibody is stable for at least two weeks at 2-4  $^{\circ}$  C.

MIDIO 44 其 数

PubMed: PubMed

产品介绍 background:

Glyceraldehyde 3 phosphate dehydrogenase (GAPDH) is well known as one of the key enzymes involved in

glycolysis. As well as functioning as a glycolytic enzyme in cytoplasm, recent evidence suggests that mammalian

GAPDH is also involved in a great number of intracellular proceses such as membrane fusion, microtubule

bundling, phosphotransferase activity, nuclear RNA export, DNA replication, and DNA repair. During the last

decade a lot of data appeared concerning the role of GAPDH in different pathologies including prostate cancer

progression, programmed neuronal cell death, age related neuronal diseases, such as Alzheimer's and

Huntington's disease. GAPDH is expressed in all cells. It is constitutively expressed in almost all tissues at high

levels. There are however some physiological factors such as hypoxia and diabetes that increase GAPDH

expression in certain cell types. GAPDH molecule is composed of four 38kDa subunits.

**Function:** 

Has both glyceraldehyde-3-phosphate dehydrogenase and nitrosylase activities, thereby playing a role in

glycolysis and nuclear functions, respectively. Participates in nuclear events including transcription, RNA

transport, DNA replication and apoptosis. Nuclear functions are probably due to the nitrosylase activity that

mediates cysteine S-nitrosylation of nuclear target proteins such as SIRT1, HDAC2 and PRKDC. Glyceraldehyde-3-

phosphate dehydrogenase is a key enzyme in glycolysis that catalyzes the first step of the pathway by converting

D-glyceraldehyde 3-phosphate (G3P) into 3-phospho-D-glyceroyl phosphate.

Subunit:

Homotetramer. Interacts with TPPP; the interaction is direct. Interacts (when S-nitrosylated) with SIAH1; leading

to nuclear translocation. Interacts with RILPL1/GOSPEL, leading to prevent the interaction between GAPDH and

SIAH1 and prevent nuclear translocation. Interacts with EIF1AD, USP25, PRKCI and WARS.

**Subcellular Location:** 

Cytoplasm, cytosol. Nucleus. Cytoplasm, perinuclear region. Membrane. Note=Translocates to the nucleus



following S-nitrosylation and interaction with SIAH1, which contains a nuclear localization signal. Postnuclear and Perinuclear regions.

Post-translational modifications:
S-nitrosylation of Cys-152 leads to interaction with SIAH1, followed by translocation to the nucleus.
ISGylated (Probable).
Sulfhydration at Cys-152 increases catalytic activity.
Similarity:
Belongs to the glyceraldehyde-3-phosphate dehydrogenase family.
SWISS:
P04406
Gene ID:
100009074
Hairana 25040 Zahasfiah
Unigene: 35640 Zebrafish
Important Note:
This product as supplied is intended for research use only, not for use in human, therapeutic or diagnostic
applications.



甘油醛-3-磷酸脱氢酶(Glyceraldehyde 3 phosphate dehydrogenase,GAPDH)是糖酵解(glycolysis)过程中的关键酶。除了在胞质中作为糖酵解的酶以外,有证据表明哺乳动物细胞中的 GAPDH 参与了多种胞内生化过程,包括膜融合(membrane fusion)、微管成束(microtubule bundling)、磷酸转移酶(phosphotransferase)激活、核内 RNA 出核、DNA 复制与 DNA 修复。一些生理因素,诸如低氧(hypoxia)和尿糖(diabetes),可以增加 GAPDH 在特定细胞中的表达。GAPDH 存在于几乎所有的组织中,以高水平持续表达。 GAPDH(甘油醛-3-磷酸脱氢酶)是参与糖酵解的一种关键酶,由 4 个 30 — 40kDa 的亚基组成。

## 产品图片:

