

## 3 号染色体开放阅读框 77 抗体

产品货号： mlR17322

英文名称： C3orf77/TOPAZ1

中文名称： 3 号染色体开放阅读框 77 抗体

别名： C3orf77; Chromosome 3 open reading frame 77; FLJ36157; Testis and ovary specific PAZ domain containing 1; testis and ovary-specific PAZ domain gene 1; Testis- and ovary-specific PAZ domain-containing protein 1; TOPAZ 1; TOPAZ1; TOPZ1\_HUMAN.

研究领域： 细胞生物 免疫学

抗体来源： Rabbit

克隆类型： Polyclonal

交叉反应： Human, Mouse, Rat, Cynomolgus Monkey, Rhesus monkey, Gorilla

产品应用： ELISA=1:500-1000 IHC-P=1:400-800 IHC-F=1:400-800 ICC=1:100-500 IF=1:100-500 (石蜡切片需

做抗原修复)

not yet tested in other applications.

optimal dilutions/concentrations should be determined by the end user.

分子量：191kDa

细胞定位：细胞浆

性状：Lyophilized or Liquid

浓度：1mg/ml

免疫原：KLH conjugated synthetic peptide derived from human C3orf77/TOPAZ1:651-750/1692

亚型：IgG

纯化方法：affinity purified by Protein A

储存液：0.01M TBS(pH7.4) with 1% BSA, 0.03% Proclin300 and 50% Glycerol.

保存条件：Store at -20 ° C for one year. Avoid repeated freeze/thaw cycles. The lyophilized antibody is stable at room temperature for at least one month and for greater than a year when kept at -20 ° C. When reconstituted in sterile pH 7.4 0.01M PBS or diluent of antibody the antibody is stable for at least two weeks at 2-4 ° C.



**PubMed :** PubMed