

Bcl-2 样蛋白 2 抗体

产品货号： mlR5903

英文名称： Bcl2L2

中文名称： Bcl-2 样蛋白 2 抗体

别名： Apoptosis regulator BCL W; Apoptosis regulator Bcl-W; B2CL2_HUMAN; BCL 2 Like 2; Bcl 2 like 2 protein; Bcl 2L2; BCL W; Bcl-2-like protein 2; Bcl2 L2; BCL2 like 2; BCL2 like 2 protein; Bcl2-L-2; Bcl2l2; BCLW.

研究领域： 肿瘤 细胞生物 免疫学 细胞凋亡

抗体来源： Rabbit

克隆类型： Polyclonal

交叉反应： Human, Mouse, Rat, Dog, Pig, Horse, Rabbit,

产品应用： WB=1:500-2000 ELISA=1:500-1000 IHC-P=1:400-800 IHC-F=1:400-800 IF=1:100-500 （石蜡切片需做抗原修复）

not yet tested in other applications.

optimal dilutions/concentrations should be determined by the end user.

分子量：21kDa

细胞定位：细胞浆 细胞膜

性状：Lyophilized or Liquid

浓度：1mg/ml

免疫原：KLH conjugated synthetic peptide derived from human Bcl2L2:1-100/193

亚型：IgG

纯化方法：affinity purified by Protein A

储存液：0.01M TBS(pH7.4) with 1% BSA, 0.03% Proclin300 and 50% Glycerol.

保存条件：Store at -20 °C for one year. Avoid repeated freeze/thaw cycles. The lyophilized antibody is stable at room temperature for at least one month and for greater than a year when kept at -20°C. When reconstituted in sterile pH 7.4 0.01M PBS or diluent of antibody the antibody is stable for at least two weeks at 2-4 °C.

PubMed：PubMed

产品介绍：Promotes cell survival. Blocks dexamethasone-induced apoptosis. Mediates survival of postmitotic Sertoli cells by suppressing death-promoting activity of BAX.

Subcellular Location:

Mitochondrion membrane; Peripheral membrane protein.

Tissue Specificity:

Expressed (at protein level) in a wide range of tissues with highest levels in brain, spinal cord, testis, pancreas, heart, spleen and mammary glands. Moderate levels found in thymus, ovary and small intestine. Not detected in salivary gland, muscle or liver. Also expressed in cell lines of myeloid, fibroblast and epithelial origin. Not detected in most lymphoid cell lines.

Similarity:

Belongs to the Bcl-2 family.

SWISS:

Q92843

Gene ID:

599

Important Note:

This product as supplied is intended for research use only, not for use in human, therapeutic or diagnostic applications.

产品图片

