

吡哆醛依赖脱羧酶结构域 PDXDC1 抗体

产品货号： mlR12606

英文名称： PDXDC1

中文名称： 吡哆醛依赖脱羧酶结构域 PDXDC1 抗体

别名： 2210010A19Rik; AA415817; KIAA0251; Kiaa0251 hp; LP8165; MGC189517; MGC8182; PDXD1_HUMAN; PDXDC 1; Pdxdc1; Pyridoxal dependent decarboxylase domain containing 1; Pyridoxal-dependent decarboxylase domain-containing protein 1; rCG_49895; RGD1562597.

研究领域： 肿瘤 细胞生物 信号转导 细胞骨架

抗体来源： Rabbit

克隆类型： Polyclonal

交叉反应： Human, Mouse, Rat, Dog, Cow, Horse, Rabbit, Sheep,

产品应用： ELISA=1:500-1000 IHC-P=1:400-800 IHC-F=1:400-800 ICC=1:100-500 IF=1:100-500 （石蜡切片需做抗原修复）

not yet tested in other applications.

optimal dilutions/concentrations should be determined by the end user.

分子量： 87kDa

细胞定位： 细胞浆

性状： Lyophilized or Liquid

浓度： 1mg/ml

免 疫 原： KLH conjugated synthetic peptide derived from human PDXDC1:651-750/788

亚 型： IgG

纯化方法： affinity purified by Protein A

储 存 液： 0.01M TBS(pH7.4) with 1% BSA, 0.03% Proclin300 and 50% Glycerol.

保存条件： Store at -20 ° C for one year. Avoid repeated freeze/thaw cycles. The lyophilized antibody is stable at room temperature for at least one month and for greater than a year when kept at -20° C. When reconstituted in sterile pH 7.4 0.01M PBS or diluent of antibody the antibody is stable for at least two weeks at 2-4 ° C.

PubMed： PubMed