

AU1 tag 标签抗体

产品货号： mlR2114

英文名称： AU1 tag

中文名称： AU1 tag 标签抗体

别名： AU1 epitope tag; DTYRYI epitope tag; DTYRYI tag.

产品类型： 标签抗体

研究领域： 免疫学

抗体来源： Rabbit

克隆类型： Polyclonal

交叉反应： AU1 tagged fusion protein

产品应用： WB=1:500-2000 ELISA=1:500-1000 ICC=1:100-500

not yet tested in other applications.

optimal dilutions/concentrations should be determined by the end user.

性状： Lyophilized or Liquid

浓度： 1mg/ml

免疫原： KLH conjugated DTYRYI epitope tag:

亚型： IgG

纯化方法： affinity purified by Protein A

储存液： 0.01M TBS(pH7.4) with 1% BSA, 0.03% Proclin300 and 50% Glycerol.



保存条件 : Store at -20 ° C for one year. Avoid repeated freeze/thaw cycles. The lyophilized antibody is stable at room temperature for at least one month and for greater than a year when kept at -20° C. When reconstituted in sterile pH 7.4 0.01M PBS or diluent of antibody the antibody is stable for at least two weeks at 2-4 ° C.

PubMed : PubMed

产品介绍 : AU5 is a commonly used epitope tag engineered onto the N- or C- terminus of a protein of interest so that the tagged protein can be analyzed and visualized. Because of the small size of the epitope, it is unlikely to alter the activity of the cloned sequence.

SWISS:

N/A

Gene ID:

N/A

Important Note:

This product as supplied is intended for research use only, not for use in human, therapeutic or diagnostic applications.