

磷酸化蛋白激酶底物相关蛋白抗体

产品货号: mlR16681

英文名称: phospho-HEF1 (Ser369)

中文名称: 磷酸化蛋白激酶底物相关蛋白抗体

别 名: HEF1 (phospho S369); p-HEF1 (phospho S369); Cas like docking; Cas scaffolding protein family member 2; CAS-L; CAS2; CasL; CASL_HUMAN; CASS2; Crk associated substrate related; Crk associated substrate related protein; CRK-associated substrate-related protein; dJ49G10.2 (Enhancer of Filamentation 1 (HEF1)); dJ761I2.1; Enhancer of filamentation 1; Enhancer of filamentation 1 p55; HEF 1; HEF1

产品类型: 磷酸化抗体

研究领域: 肿瘤 细胞生物 信号转导 细胞凋亡 转录调节因子 表观遗传学

抗体来源: Rabbit

克隆类型: Polyclonal



交叉反应: Human, Cow, Horse, Sheep,

产品应用: WB=1:500-2000 ELISA=1:500-1000 IHC-P=1:400-800 IHC-F=1:400-800 ICC=1:100-500 IF=1:100-500 (石蜡切片需做抗原修复)

not yet tested in other applications.

optimal dilutions/concentrations should be determined by the end user.

分子量: 93kDa

细胞定位: 细胞浆

性 状: Lyophilized or Liquid

浓 度: 1mg/ml

免 疫 原 : KLH conjugated synthesised phosphopeptide derived from human HEF1 around the phosphorylation site of Ser369:RL(p-S)FS

亚 型: IgG

纯化方法: affinity purified by Protein A

储存液: 0.01M TBS(pH7.4) with 1% BSA, 0.03% Proclin300 and 50% Glycerol.

mbio 4数数数

保存条件: Store at -20 °C for one year. Avoid repeated freeze/thaw cycles. The lyophilized antibody is

stable at room temperature for at least one month and for greater than a year when kept at -20° C. When

reconstituted in sterile pH 7.4 0.01M PBS or diluent of antibody the antibody is stable for at least two weeks at 2-

4 ° C.

PubMed: PubMed

产品介绍 background:

The protein encoded by this gene is a member of the CRK-associated substrates family. Members of this family

are adhesion docking molecules that mediate protein-protein interactions for signal transduction pathways. This

protein is a focal adhesion protein that acts as a scaffold to regulate signaling complexes important in cell

attachment, migration and invasion as well as apoptosis and the cell cycle. This protein has also been reported to

have a role in cancer metastasis. Alternative splicing results in multiple transcript variants. [provided by RefSeq,

Aug 2012]

Function:

Docking protein which plays a central coordinating role for tyrosine-kinase-based signaling related to cell

adhesion. May function in transmitting growth control signals between focal adhesions at the cell periphery and

the mitotic spindle in response to adhesion or growth factor signals initiating cell proliferation. May play an

important role in integrin beta-1 or B cell antigen receptor (BCR) mediated signaling in B- and T-cells. Integrin

beta-1 stimulation leads to recruitment of various proteins including CRK, NCK and SHPTP2 to the tyrosine

phosphorylated form.

Subcellular Location:

Cytoplasm; cytoskeleton; spindle and Cytoplasm; cell cortex. Nucleus. Golgi apparatus. Cell projection;

lamellipodium. Cytoplasm. Cell junction; focal adhesion. Localizes to both the cell nucleus and the cell periphery

and is differently localized in fibroblasts and epithelial cells. In fibroblasts is predominantly nuclear and in some

cells is present in the Golgi apparatus. In epithelial cells localized predominantly in the cell periphery with

particular concentration in lamellipodia but is also found in the nucleus. Isoforms p105 and p115 are

predominantly cytoplasmic and associate with focal adhesions while p55 associates with mitotic spindle.



Tissue Specificity:

Widely expressed. Higher levels detected in kidney, lung, and placenta. Also detected in T-cells, B-cells and

diverse cell lines. The protein has been detected in lymphocytes, in diverse cell lines, and in lung tissues.

Post-translational modifications:

Cell cycle-regulated processing produces four isoforms: p115, p105, p65, and p55. Isoform p115 arises from p105

phosphorylation and appears later in the cell cycle. Isoform p55 arises from p105 as a result of cleavage at a

caspase cleavage-related site and it appears specifically at mitosis. The p65 isoform is poorly detected. Focal

adhesion kinase 1 phosphorylates the protein at the YDYVHL motif (conserved among all cas proteins). The SRC

family kinases (FYN, SRC, LCK and CRK) are recruited to the phosphorylated sites and can phosphorylate other

tyrosine residues. Ligation of either integrin beta-1 or B-cell antigen receptor on tonsillar B-cells and B-cell lines

 $promotes\ tyrosine\ phosphorylation\ and\ both\ integrin\ and\ BCR-mediated\ tyrosine\ phosphorylation\ requires\ an$

intact actin network. In fibroblasts transformation with oncogene v-ABL results in an increase in tyrosine phosphorylation. Transiently phosphorylated following CD3 cross-linking and this phosphorylated form binds to

CRK and C3G. A mutant lacking the SH3 domain is phosphorylated upon CD3 cross-linking but not upon integrin

beta-1 cross-linking. Tyrosine phosphorylation occurs upon stimulation of the G-protein coupled C1a calcitonin

receptor in rabbit. Calcitonin-stimulated tyrosine phosphorylation is mediated by calcium- and protein kinase C-

dependent mechanisms and requires the integrity of the actin cytoskeleton.

Similarity:

Belongs to the CAS family.

Contains 1 SH3 domain.

SWISS:

Q14511

Gene ID:



4739

Important Note:

This product as supplied is intended for research use only, not for use in human, therapeutic or diagnostic applications.