

磷酸化蛋白酪氨酸磷酸酶 1C 抗体

产品货号: mlR10167

英文名称: phospho-PTPN6 (Tyr564)

中文名称: 磷酸化蛋白酪氨酸磷酸酶 1C 抗体

别 名: SHP-1(Phospho-Tyr564); PTPN6(phospho Y564); SHP1(phospho Y564); SHP1; 70 kda SHP 1L protein; 70 kda SHP1L protein; 70Z-SHP; EC 3.1.3.48; HCP; HCPH; Hematopoietic cell phosphatase; Hematopoietic cell protein tyrosine phosphatase; Hematopoietic cell protein-tyrosine phosphatase; HPTP 1C; HPTP1C; MGC124580; Protein tyrosine phosphatase 1C; Protein tyrosine phosphatase non receptor type 6; Protein tyrosine phosphatase SHP 1; Protein tyrosine phosphatase SHP1; Protein-tyrosine phosphatase 1C; Protein-tyrosine phosphatase SHP-1; PTN6_HUMAN; PTP 1C; PTP-1C; PTPN 6; PTPN6; SH PTP 1; SH PTP1; SH-PTP1; SHP 1; SHP1L; SHPTP 1; SHPTP1; Tyrosine protein phosphatase non receptor type 6; Tyrosine-protein phosphatase non-receptor type 6.

产品类型: 磷酸化抗体

抗体来源: Rabbit

克隆类型: Polyclonal

交叉反应: Human, Mouse, Rat, Dog, Pig, Horse, Sheep, Guinea Pig,

产品应用 : WB=1:500-2000 ELISA=1:500-1000 IHC-P=1:400-800 IHC-F=1:400-800 ICC=1:100-500 IF=1:100-500

(石蜡切片需做抗原修复)

not yet tested in other applications.

optimal dilutions/concentrations should be determined by the end user.

分子量: 65kDa

细胞定位: 细胞核 细胞浆



性 状: Lyophilized or Liquid

浓 度: 1mg/ml

兔 疫 原 : KLH conjugated Synthesised phosphopeptide derived from human PTPN6 around the

phosphorylation site of Tyr564:DV(p-Y)EN

亚型: IgG

纯化方法: affinity purified by Protein A

储存液: 0.01M TBS(pH7.4) with 1% BSA, 0.03% Proclin300 and 50% Glycerol.

保存条件: Store at -20 °C for one year. Avoid repeated freeze/thaw cycles. The lyophilized antibody is stable

at room temperature for at least one month and for greater than a year when kept at -20°C. When reconstituted

in sterile pH 7.4 0.01M PBS or diluent of antibody the antibody is stable for at least two weeks at 2-4 °C.

PubMed: PubMed

产品介绍: The protein encoded by this gene is a member of the protein tyrosine phosphatase (PTP) family.

PTPs are known to be signaling molecules that regulate a variety of cellular processes including cell growth,

differentiation, mitotic cycle, and oncogenic transformation. N-terminal part of this PTP contains two tandem Src

homolog (SH2) domains, which act as protein phospho-tyrosine binding domains, and mediate the interaction of

this PTP with its substrates. This PTP is expressed primarily in hematopoietic cells, and functions as an important

regulator of multiple signaling pathways in hematopoietic cells. This PTP has been shown to interact with, and

dephosphorylate a wide spectrum of phospho-proteins involved in hematopoietic cell signaling. Multiple

alternatively spliced variants of this gene, which encode distinct isoforms, have been reported. [provided by

RefSeq, Jul 2008]

Function:

Modulates signaling by tyrosine phosphorylated cell surface receptors such as KIT and the EGF receptor/EGFR.

The SH2 regions may interact with other cellular components to modulate its own phosphatase activity against

interacting substrates. Together with MTUS1, induces UBE2V2 expression upon angiotensin II stimulation. Plays a

key role in hematopoiesis.



Subunit:

Monomer. Interacts with MTUS1. Interacts with MILR1 (tyrosine-phosphorylated). Interacts with KIT. Binds PTPNS1, LILRB1 and LILRB2. Interacts with FCRL2, FCRL3, FCRL4, CD300LF, CDK2 and CD84. Interacts with KIR2DL1; the interaction is enhanced by ARRB2. Interacts (via SH2 1 domain) with ROS1; the interaction is direct and promotes ROS1 dephosphorylation. Interacts with EGFR; inhibits EGFR-dependent activation of MAPK/ERK. Interacts with LYN. Interacts with the tyrosine phosphorylated form of PDPK1.

Subcellular Location:

Cytoplasm. Nucleus. Note=In neurons, translocates into the nucleus after treatment with angiotensin II. Shuttles between the cytoplasm and nucleus via its association with PDPK1.

Tissue Specificity:

Isoform 1 is expressed in hematopoietic cells. Isoform 2 is expressed in non-hematopoietic cells.

Post-translational modifications:

Phosphorylated on tyrosine residues. Binding of KITLG/SCF to KIT increases tyrosine phosphorylation. Phosphorylation at Tyr-564 enhances phosphatase activity.

Similarity:

Belongs to the protein-tyrosine phosphatase family. Non-receptor class 2 subfamily.

Contains 2 SH2 domains.

Contains 1 tyrosine-protein phosphatase domain.

SWISS:

P29350



_		
ϵ	20	IID.
GE.	пе	II.

5777

Important Note:

This product as supplied is intended for research use only, not for use in human, therapeutic or diagnostic applications.

产品图片

