

细胞分裂周期 2 样蛋白激酶 6 抗体

产品货号： mlR7070

英文名称： CDC2L6

中文名称： 细胞分裂周期 2 样蛋白激酶 6 抗体

别名： CDC2 related protein kinase 6; CDC2-related protein kinase 6; CDK11; Cdk19; CDK19_HUMAN; CDK8 like cyclin dependent kinase; Cell division cycle 2 like 6 (CDK8 like); Cell division cycle 2 like protein kinase 6; Cell division cycle 2-like protein kinase 6; Cell division protein kinase 19; Cyclin dependent kinase (CDC2 like) 11; Cyclin dependent kinase 11; Cyclin-dependent kinase 11; Cyclin-dependent kinase 19; Death preventing kinase; Death-preventing kinase.

研究领域： 细胞生物 信号转导 细胞周期蛋白 激酶和磷酸酶 表观遗传学

抗体来源： Rabbit

克隆类型： Polyclonal

交叉反应： Human, Mouse, Rat, Cow, Horse, Rabbit, Monkey, Chimpanzee,

产品应用： ELISA=1:500-1000 IHC-P=1:400-800 IHC-F=1:400-800 IF=1:100-500 (石蜡切片需做抗原修复)

not yet tested in other applications.

optimal dilutions/concentrations should be determined by the end user.

分子量：57kDa

细胞定位：细胞核 细胞浆

性状：Lyophilized or Liquid

浓度：1mg/ml

免疫原：KLH conjugated synthetic peptide derived from human CDC2L6:21-120/502

亚型：IgG

纯化方法：affinity purified by Protein A

储存液：0.01M TBS(pH7.4) with 1% BSA, 0.03% Proclin300 and 50% Glycerol.

保存条件：Store at -20 °C for one year. Avoid repeated freeze/thaw cycles. The lyophilized antibody is stable at room temperature for at least one month and for greater than a year when kept at -20°C. When reconstituted in sterile pH 7.4 0.01M PBS or diluent of antibody the antibody is stable for at least two weeks at 2-4 °C.

PubMed：PubMed

产品介绍： This gene encodes a protein that is one of the components of the Mediator coactivator complex. The Mediator complex is a multiprotein complex required for transcriptional activation by DNA binding transcription factors of genes transcribed by RNA polymerase II. The protein encoded by this gene is similar to cyclin-dependent kinase 8 which can also be a component of the Mediator complex.

Similarity:

Belongs to the protein kinase superfamily. CMGC Ser/Thr protein kinase family. CDC2/CDKX subfamily.

Contains 1 protein kinase domain.

SWISS:

Q9BWU1

Gene ID:

23097

Important Note:

This product as supplied is intended for research use only, not for use in human, therapeutic or diagnostic applications.

产品图片

